



PATHOLOGY WORKSTATION BAIP-301

PATHOLOGY WORKSTATION BAIP-301

Workstations are used during dissections, autopsies, necropsies, and other pathological procedures. It includes grossing/trimming stations and dissection tables. It may also include a built-in sink and ventilation system.

Used in Hospital pathology department, Pathology laboratory.

Also known as Pathological sampling bench.

BAIP-301 PATHOLOGY WORKSTATION



Odor-remove function

Single water sink, two-person operations

Exhaust air system: High-power suction device, low noise

Optional water heater is available, realizing hot & cold water adjustment

The base cabinet is inward-concave type, convenient for sitting operation

Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes

Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning

SPECIFICATIONS

| Model | BAIP-301 |
|---------------------------|-------------------------------------|
| Work Surface Height | 850 mm |
| Exhaust Airflow Volume | 500~600m ³ /h |
| Main Body Material | 304 stainless steel |
| Work Table Material | 304 stainless steel |
| Package Size (WxDxH) | 1160x890x2170 mm |
| Gross Weight | 250 kg |
| Working Area Size (WxDxH) | 895x430x900 mm |
| External Size (WxDxH) | 1000x750x1950 mm |
| Consumption | 250 W |
| Power Supply | AC220V±10%, 50/60Hz; 110V±10%, 60Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com