



## SLIDE-DRYING BOARDS BCBT-304

## SLIDE-DRYING BOARDS BCBT-304

A metal mould is specifically used to create the slots. The stamping and shaping processes are linked to create the labelling slots.

Used in Pathological laboratories, Research institutions, Universities, Cytology and Histology applications.

## BCBT-304 SLIDE-DRYING BOARDS



The slots are specifically manufactured using a metal mold.

The labelling slots are manufactured through an integrated stamping and forming operation.

## SPECIFICATIONS

Model	BCBT-304
Capacity	96 microscope slide-drying boards (24 in each section)
Base Material	SPCC1.2 cold-roll steel sheets
Main body and Drawers Material	SPCC 0.8 cold-roll steel sheets
Dimensions	450x315x1210 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)