





# LOW TEMPERATURE CIRCULATOR BCLT-2202





### LOW TEMPERATURE CIRCULATOR BCLT-2202

## **BCLT-2202 LOW TEMPERATURE CIRCULATOR**



P.I.D temperature controller provides accurate and reliable temperature control Large LCD display screen and interface provides for user-friendly operation Preset On/Off function

Independent circulating pump switch for easy starting/stopping the cycle Easy to set Adjustable Timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel No angle in bath chamber, easy to clean. Uniquely designed heating method can up to 80  $^\circ\text{C}$ 

Programming setting function with 7 periods and 9 steps for each period, which means. There are 63 programmable steps in total

## **SPECIFICATIONS**

Model	BCLT-2202
Capacity	8 L
Temperature range	-20 - 80°C
Temperature accuracy	±0.5°C
Cooling capacity at 10°C	600 W
Refrigerant	R404A
Security features	Delay, leakage, overcurrent, overvoltage
Total power	1350 W
Power requirements	AC 220 V / 50 HZ
Pump power	30
Pump flow max.	6 L / min
Maximum head	9 m
Inlet/Outlet pipe diameter	Ф10 mm
Weight	41
Liquid storage	5.2
Noise level	≤45
Dimension	310Wx445Dx570H mm



### Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada Email: info@biolabscientific.com | Website: www.biolabscientific.com