


## PRODUCT CATALOG



# **ULTRA LOW RETENTION TIP BPIC-408**



www.biolabscientific.com

### ULTRA LOW RETENTION TIP BPIC-408

Pipette tips are essential when using a micropipette. These get attached to the shaft of the pipette and are usually ejected after use. It is an essential lab item for any test or experiment. They are available as racked or loose bulk tips, with or without a filter, sterile or non-sterile.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Pipet Tips, Laboratory Pipette Tips, Laboratory Pipet Tips.

#### **BPIC-408 ULTRA LOW RETENTION TIP**



Made from medical grade Polypropylene Slippy surface with excellent design DNase-free Rnase-free Non-pyrogenic Filter and non-filter for choice Bag and rack are available to pack

#### SPECIFICATIONS

Model	BPIC-408
Volume	200 µl
Package	1000 pcs/bag
Unit	10



Biolab Scientific Ltd. 3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada Email: info@biolabscientific.com | Website: www.biolabscientific.com