

# PRODUCT CATALOG



## **BENCHTOP SHAKING INCUBATOR BSBT-2103**





#### **BENCHTOP SHAKING INCUBATOR BSBT-2103**

#### **BSBT-2103 BENCHTOP SHAKING INCUBATOR**



Adopt 4.3 inches touch screen, intelligent touch screen control, parameters real time display, operation easily

Temp., rotation speed and time etc parameters fast setting

Silent working.

Colorful touch screen display, multi data display at one screen, easily to observe and operate

Platform is stainless steel, easy to detach and clean

Various kinds of clamps and racks for selection, easily change racks and clamps improve work efficient

User's set parameters can be stored automatically when power off suddenly, and it resume last program settings when power on

PID microcomputer control temp. and shake frequency with timer, it starts shake gently to prevent liquid coming out from the vessels and stop the machine working Brushless DC motor, large start torque, wide speed adjustment, free maintenance

Three-eccentric axis balance drive ensures all samples in the shaker platform

movements at the same rotation frequency. Durable structure guarantees the shaker incubator can work properly no matter full load or high speed working

Multi-segment temp., speed, time program at the same time. It has normal mode or programmable mode for running. Preset parameter and work data can be shown at the same time. It simplified complicated culture requirements and realize auto control and run

Equip with USB as options, record temp. parameter changes

Equip with RS232 data interface, through software to remote control machine performance

### **SPECIFICATIONS**

BSBT-2103
40-300 rpm
26 mm
1 - 99h59min
250mlx63pcs or 500mlx35pcs or 1000mlx35pcs or 2000mlx16pcs (single shelf)
920x640 mm
1 pcs
1115Wx835Dx460H mm
250 W
AC 220 V 50 HZ
1000mlx35pcs (single shelf)



#### Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada Email: info@biolabscientific.com | Website: www.biolabscientific.com