



## ROTATOR BTUR-202

## ROTATOR BTUR-202

The versatile, end-over-end tube roller and Rotators are designed for mixing samples in both a vertical and horizontal plane. They can be placed in an environment temperature ranges from 5 to 40° and the rotation speed is adjustable, from 0 to 99 rpm. Tube Rollers, Rotators and Rocker are designed to deliver salient performance, accuracy, and reproducible results for all your application needs.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## BTUR-202 ROTATOR



The rotator comes with a disk type design  
Digital control speed and allowing procedures to be accurate  
LCD displays timer  
End-over-end or rolling action, gives gentle but effective mixing  
Adjustable speed range of 10 to 70rpm  
Adjustable mixing angle between 0 and 90 °  
Used with a wide variety of accessories  
Continuous and timed quiet operation

## SPECIFICATIONS

Model	BTUR-202
Speed Range	10-70 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	0.8
Tilt Angle (°)	0-90
Operation	Continuous/ timed operation
Timer	Yes
Timing Range	1-1199 min
Dimension	510Dx190Wx260H mm
Motor Type	DC motor
Protection Class	IP21
Display	LCD
Weight	8 kg
Power	40 W
Power Supply	100-240V, 50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3300807006	Disk accessory	For 1.5ml x 60 centrifuge tubes holder
3300807007	Disk accessory	For 15ml x 16 centrifuge tubes holder
3300807008	Disk accessory	For 50ml x 8 centrifuge tubes holder

3300807009

Disk support rods

4 pcs; use with BX-103920; BX-103921; BX-103922



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)