



Product Image Coming Soon

MEDICAL AUTOCLAVE

MEDICAL AUTOCLAVE

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with door locking system, these autoclaves are ideal for small spaces.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

400 MEDICAL CLASS N AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	BAMN-401	BAMN-403	BAMN-404
Capacity	3 L	18 L	23 L
Sterilization Pressure	0.22 Mp	-	-
Sterilization Temperature	134 °C	-	-
Sterilization Time	4 min	-	-
Sterilization Cycle time	11 min	-	-
Required water	Distilled water or purified water	-	-
Product Size	300X190X350 mm	450X550X350 mm	450x550x350 mm
Outer Dimension	370x260x460 mm	-	-
Gross Weight	7.5 Kg	30 Kg	33 Kg
Power	0.5 KW	800 W	
Voltage	AC220/110V+10%	AC220/110V+10%, 50/60 Hz+ 1Hz	
Chamber Dimension	-	φ249x350 mm	φ249x450 mm
Package Size	-	670x550x500 mm	
Net Weight	-	23 Kg	26 Kg



Product Image Coming Soon

BAMN-401



Product Image Coming Soon

BAMN-403



Product Image Coming Soon

BAMN-404



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com