



LABORATORY VERTICAL AUTOCLAVE

LABORATORY VERTICAL AUTOCLAVE

We offer a variety of vertical laboratory autoclaves for general purpose sterilization applications. Customers can choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Food industry, Chemical industry, Research, Laboratory, Pharmaceutical.

Also known as Vertical Pressure Sterilizer, Vertical Pressure Autoclave, Vertical Sterilizer, Laboratory Vertical Autoclave,.

100 VERTICAL AUTOCLAVE



- Stainless steel structure
- Hand wheel type / Screw type
- Door safety lock system
- Digital Display
- Automatic exhaust system for cold air and pressure
- Over-temperature and over-pressure auto-protection
- Alarm function for completion of process and low water level
- Self-inflating type seal
- Optional: Drying system

SPECIFICATIONS

Model	BAVT-101	BAVT-102	BAVT-103	BAVT-104
Capacity	35L	50 L	75 L	100 L
Temperature Range	105~134°C			
Heat Average	$\pm 1^{\circ}\text{C}$			
Working Pressure	0.22 MPa			
Maximum Working Pressure	0.23 Mpa			
Working Temperature	134 °C			
Timer Range	0~90 min or 0~99hour59min			
Overall Dimension (LxWxH)	470x450x1020 mm	520x490x1150 mm	550x530x1180 mm	580x560x1290 mm
Packing Size	540x540x1140 mm	580x580x1250 mm	610x610x12580 mm	640x640x1360 mm
Weight	81 / 63 kg	94 / 74 kg	113 / 93 kg	125 / 105 kg
Power Supply	1x2.5KW/AC220V 50Hz	2x1.5KW/AC220V 50Hz	3x1.5KW/AC220V 50Hz	



BAVT-101



BAVT-102



BAVT-103



BAVT-104

300 VERTICAL AUTOCLAVE



Microcomputer control system
 Unique anti scald design chamber for corrosion resistance
 Stainless steel chamber
 Five sterilization modes
 User friendly control panel
 Multi layered pressure protection
 Fully automatic internal discharge control embedded with a steam trapping bottle
 Automatic checking of locking condition of cover
 Automatic sterilization according to preset time
 Safety device: Self-induction pressure interlocking device, over temperature protection system, dry scorch protection system, overpressure protection, safety valve, over current and short circuit protection, cover inspection system, leakage protection device, anti-scald safety protection and automatic troubleshooting system

SPECIFICATIONS

Model	BAVT-302-A	BAVT-302-B	BAVT-303-A
Capacity	54 L		80 L
Function	Fully automatic, Drying function	Fully automatic	Fully automatic, Drying function
Temperature Range	5°C-138°C		
Sterilization Temperature	105-135°C		
Warming Temperature	45-60°C		
Warming Time	1-9999 min		
Melting Temperature	60-100		
Melting Time	1-300 min		
Humidity Range	10%-85% RH		
Working Pressure	0-0.6 MPa		
Chamber Dimension	320x730 mm		400x700 mm
Overall Dimension (LxWxH)	460x542x1070 mm		620x680x1040 mm
Weight	80/120 kg		140/180 kg
Power	2900 W		4600 W
Power Supply	220 V, 60Hz		

Model	BAVT-303-B	BAVT-304-A	BAVT-304-B
Capacity	80 L	100 L	
Function	Fully automatic	Fully automatic, Drying function	Fully automatic
Temperature Range	5°C-138°C		
Sterilization Temperature	105-135°C		
Warming Temperature	45-60°C		
Warming Time	1-9999 min		
Melting Temperature	60-100		
Melting Time	1-300 min		
Humidity Range	10%-85% RH		
Working Pressure	0-0.6 MPa		
Chamber Dimension	400x700 mm	400x800 mm	
Overall Dimension (LxWxH)	620x680x1040 mm	620x680x1140 mm	

Weight	140/180 kg	145/185 kg
Power	4600 W	
Power Supply	220 V, 60Hz	

ACCESSORIES

Accessory Code	Name
1702806006	Stainless steel baskets
1702806007	Drainage hose
1702806008	Exhaust bottle

OPTIONAL ACCESSORIES

Accessory Code	Name
1702806009	Printer
1702806010	Printing set
1702806011	Load thermometer
1702806012	Cooling fan
1702806013	Temperature controller
1702806014	Temperature sensor
1702806015	Pressure Sensor
1702806016	Heater
1702806017	Sealing ring

Accessory Code	Name
1702806018	Pressure gauge
1702806019	Stainless steel pail
1702806020	Relay
1702806021	Main board
1702806022	Display board
1702806023	Jiggle switch
1702806024	Handle
1702806025	Pressure switch
1702806026	Solenoid Valve



BAVT-302-A



BAVT-302-B



BAVT-303-A



BAVT-303-B



BAVT-304-A



BAVT-304-B

400 VERTICAL AUTOCLAVE



- Microcomputer control system
- Stainless steel chamber
- Six sterilization modes
- Resin coated insulation of Chamber cover and bench to avoid hand burns
- User friendly control panel
- Waste sterilization mode and drying function
- Rounded corners avoid collisions and injuries
- Automatic sterilization according to preset time
- Temperature and pressure calibration interface
- Convenience of frequent water feeding
- Automatic monitoring of the discharge conditions of cold air
- Automatic opening of the vent valve to discharge steam on completion of sterilization procedure
- Fully automatic internal discharge control embedded with a steam trapping bottle
- Safety function: "Eight-column evenly distributed" interlocking device, electric double inner locks, dry scorch protection, overpressure protection, safety valve, over temperature protection, overcurrent and short circuit protection, leakage protection, cooling lock, anti-scald chamber cover and bench, automatic troubleshooting

SPECIFICATIONS

Model	BAVT-401-A	BAVT-401-B	BAVT-401-C
Capacity	60 L		
Function	Fully automatic, Automatic water feeding, Fast cooling fan, Drying function	Fully automatic, Automatic water feeding, Fast cooling fan	Fully automatic
Sterilization Temperature	105-138°C		
Sterilization Time	1-300 min		
Warming Temperature	45-60°C		
Warming Time	1-9999 min		
Melting Temperature	60-100		
Melting Time	1-300 min		
Steam Exhausting Temperature	73-104°C		
Humidity Range	10%-85% RH		
Working Pressure	0-0.6 MPa		
Maximum Working Pressure	0.27 MPa		
Timer Range	0 minute-6 days		
Chamber Dimension	400x505 mm		
Overall Dimension (LxWxH)	660x644x980 mm		
Weight	101/140 kg		
Power	2900 W		
Power Supply	220 V, 60Hz		

Model	BAVT-402-A	BAVT-402-B	BAVT-402-C
Capacity	85 L		
Function	Fully automatic, Automatic water feeding, Fast cooling fan, Drying function	Fully automatic, Automatic water feeding, Fast cooling fan	Fully automatic
Sterilization Temperature	105-138°C		

Sterilization Time	1-300 min
Warming Temperature	45-60°C
Warming Time	1-9999 min
Melting Temperature	60-100
Melting Time	1-300 min
Steam Exhausting Temperature	73-104°C
Humidity Range	10%-85% RH
Working Pressure	0-0.6 MPa
Maximum Working Pressure	0.27 MPa
Timer Range	0 minute-6 days
Chamber Dimension	400x700 mm
Overall Dimension (LxWxH)	660x644x980 mm
Weight	113/160 kg
Power	4600 W
Power Supply	220 V, 60Hz

Model	BAVT-403-A	BAVT-403-B	BAVT-403-C
Capacity	110 L		
Function	Fully automatic, Automatic water feeding, Fast cooling fan, Drying function	Fully automatic, Automatic water feeding, Fast cooling fan	Fully automatic
Sterilization Temperature	105-138°C		
Sterilization Time	1-300 min		
Warming Temperature	45-60°C		
Warming Time	1-9999 min		
Melting Temperature	60-100		
Melting Time	1-300 min		
Steam Exhausting Temperature	73-104°C		
Humidity Range	10%-85% RH		
Working Pressure	0-0.6 MPa		
Maximum Working Pressure	0.27 MPa		
Timer Range	0 minute-6 days		
Chamber Dimension	400x895 mm		
Overall Dimension (LxWxH)	660x644x1180 mm		
Weight	125/170 kg		
Power	4600 W		
Power Supply	220 V, 60Hz		


ACCESSORIES

Accessory Code	Name
1702806006	Stainless steel baskets
1702806007	Drainage hose
1702806008	Exhaust bottle

OPTIONAL ACCESSORIES

Accessory Code	Name
1702806009	Printer
1702806010	Printing set
1702806011	Load thermometer
1702806012	Cooling fan
1702806013	Temperature controller
1702806014	Temperature sensor
1702806015	Pressure Sensor
1702806016	Heater
1702806017	Sealing ring

Accessory Code	Name
1702806018	Pressure gauge
1702806019	Stainless steel pail
1702806020	Relay
1702806021	Main board
1702806022	Display board
1702806023	Jiggle switch
1702806024	Handle
1702806025	Pressure switch
1702806026	Solenoid Valve



BAVT-401-A



BAVT-401-B




BAVT-401-C



BAVT-402-A




BAVT-402-B



BAVT-402-C



BAVT-403-A



BAVT-403-B



BAVT-403-C

500 VERTICAL AUTOCLAVE



Sterilization process can be pre-set and completed automatically

Simple and easily understandable control panel to show the sterilization temperature, sterilization time and drying time at the same time, pressing separate keys to set the temperature and time makes the operation easier

System automatically monitors the discharging of the cold air to guarantee a pure steam sterilization environment and optimal sterilization effects

LED screen clearly shows the dynamic sterilization process

13 error detection procedures to monitor the normal operation of the instrument

Self-induction interlocking devices to avoid risk of misuse

Dual over-pressure protection system: upon detection of any abnormal pressure, system will relieve pressures or disconnect power supplies

Over temperature protection system: power supply will be disconnected for protection upon detection of over temperature

Dry soorch protection system: power supply will be disconnected for protection upon detection of overflow water level

Cover checking system: system can automatically check the locking of the cover, if the cover is not tightly locked, the autoclave can't be started

Electric leakage, over current and short circuit protection: system will disconnect the power supply to avoid an electricity shock or other accidents

Unique anti-scald cover and bench: avoid causing corrosions and scald

SPECIFICATIONS

Model	BAVT-502-B	BAVT-502-A	BAVT-503-B
Capacity	50 L		80 L
Function	Fully automatic	Fully automatic + Drying	Fully automatic
Sterilization chamber (Dia*mm)	Ø325x670		Ø400x709
Working environment	5°C-40°C, relative humidity 10%-85%		
Chamber material	SUS 304 stainless steel		
Sterilization Temperature	105-135°C		
Sterilization Time	1-300 min		
Range of pressure gauge	0-0.4 Mpa		
Controller	Microcomputer control system		
Mode	Solid mode, liquid mode		
Stainless steel basket (Dia*H,mm)	(Ø300x200)x3		(Ø360x280)x2
Overall Dimension (LxWxH)	460x542x1070 mm		620x680x1042 mm
Power	2600		4600
Standard Accessories	Stainless steel baskets, waste bottle, water plate, drainage hose		
Optional Accessories	Printer, printing set, stainless steel pails		
Power Supply	AC220V 50Hz		
Drying time range	-	1-300 min	-

Model	BAVT-503-A	BAVT-504-B	BAVT-504-A
Capacity	80 L	100 L	
Function	Fully automatic + Drying	Fully automatic	Fully automatic + Drying
Sterilization chamber (Dia*mm)	Ø400x709	Ø325x889	
Working environment	5°C-40°C, relative humidity 10%-85%		

Chamber material	SUS 304 stainless steel		
Sterilization Temperature	105-135℃		
Sterilization Time	1-300 min		
Range of pressure gauge	0-0.4 Mpa		
Controller	Microcomputer control system		
Mode	Solid mode, liquid mode		
Stainless steel basket (Dia*H,mm)	(Ø360x280)x2	(Ø360x220)x3	
Overall Dimension (LxWxH)	620x680x1042 mm	620x680x1222 mm	
Power	4600		
Standard Accessories	Stainless steel baskets, waste bottle, water plate, drainage hose		
Optional Accessories	Printer, printing set, stainless steel pails		
Power Supply	AC220V 50Hz		
Drying time range	1-300 min	-	1-300 min



BAVT-502-B



BAVT-502-A



BAVT-503-B



BAVT-503-A



BAVT-504-B



BAVT-504-A

600 VERTICAL AUTOCLAVE



Sterilization process can be pre-set and completed automatically

Simple and easily understandable control panel to show the sterilization temperature, sterilization time and drying time at the same time, pressing separate keys to set the temperature and time makes the operation easier

System automatically monitors the discharging of the cold air to guarantee a pure steam sterilization environment and optimal sterilization effects

LED screen clearly shows the dynamic sterilization process

13 error detection procedures to monitor the normal operation of the instrument

Self-induction interlocking devices to avoid risk of misuse

Dual over-pressure protection system: upon detection of any abnormal pressure, system will relieve pressures or disconnect power supplies

Over temperature protection system: power supply will be disconnected for protection upon detection of over temperature

Dry soorch protection system: power supply will be disconnected for protection upon detection of overflow water level

Cover checking system: system can automatically check the locking of the cover, if the cover is not tightly locked, the autoclave can't be started

Electric leakage, over current and short circuit protection: system will disconnect the power supply to avoid an electricity shock or other accidents

Unique anti-scald cover and bench: avoid causing corrosions and scald

SPECIFICATIONS

Model	BAVT-601-A	BAVT-601-B	BAVT-602-A	BAVT-602-B
Capacity	29 L		54 L	
Sterilization Temperature	105-135℃			
Sterilization Time	1-999 min			
Sterilization mode	1.Liquid 2.Liquid with warming 3.Solid 4.Instruments packed with cloth 5.Fabric 6.Rubber 7.Fast 8.Waste 9.Agar melting 10.Self defined	1.Solid 2.Agar with warming 3.Self-defined		
Melting temperature	60-100℃			
Melting time range	1-999min			
Warming Temperature	45-60℃			
Warming time range	1-9999min			
Exhaust temperature range	73-104℃(Liquid mode) 73-127℃(Solid mode)	-	73-104℃(Liquid mode) 73-127℃(Solid mode)	-
Exhaust mode	Automatic internal discharging with an embedded steam collector			
Timer range	0-6 days(Sterilization auto start up appointment)	-	0-6 days(Sterilization auto start up appointment)	-
Pressure gauge range	-0.1-0.5Mpa			
Pressure	Safety valve set pressure:2.9bar			
Controller	"Smart II" fast speed microcomputer controller			
Working environment	5-45℃,relatively humidity 10%-85%			
Chamber material	SUS304			
Chamber(Diax H)	325x381			
Rate power supply(W)	2600			
Power requirement	220V±10%,16A,50/60Hz			

Dimension(LxWxH,mm)	460x542x940		460x542x1070	
Baskets(DiaxH)	(300x200)x2		(300x200)x3	
Safety device	Seft-induction interlocking device, lid closing checking, over temperature protection system, dry scorch protection, over pressure protection, safety valve, over current and short circuit protection, electric leakage protection, scald protection, automatic failure-inspection system.			
Storage	up to 60 programs	up to 20 programs	up to 60 programs	up to 20 programs

Model	BAVT-603-A	BAVT-603-B	BAVT-604-A	BAVT-604-B
Capacity	80 L		100 L	
Sterilization Temperature	105-135℃			
Sterilization Time	1-999 min			
Sterilization mode	1.Liquid 2.Liquid with warming 3.Solid 4.Instruments packed with cloth 5.Fabric 6.Rubber 7.Fast 8.Waste 9.Agar melting 10.Self defined			
Melting temperature	60-100℃			
Melting time range	1-999min			
Warming Temperature	45-60℃			
Warming time range	1-9999min			
Exhaust temperature range	-	-	-	-
Exhaust mode	Automatic internal discharging with an embedded steam collector			
Timer range	0-6 days(Sterilization auto start up appointment)	-	0-6 days(Sterilization auto start up appointment)	-
Pressure gauge range	-0.1-0.5Mpa			
Pressure	Safety valve set pressure:2.9bar			
Controller	“Smart II” fast speed microcomputer controller			
Working environment	5-45℃,relatively humidity 10%-85%			
Chamber material	SUS304			
Chamber(Diax H)	400x709		400x809	
Rate power supply(W)	4600			
Power requirement	220V±10%,32A,50/60Hz			
Dimension(LxWxH,mm)	620x680x1100		620x680x1200	
Baskets(DiaxH)	(360x280)x2		(360x220)x3	
Safety device	Seft-induction interlocking device, lid closing checking, over temperature protection system, dry scorch protection, over pressure protection, safety valve, over current and short circuit protection, electric leakage protection, scald protection, automatic failure-inspection system.			
Storage	up to 60 programs	up to 20 programs	up to 60 programs	up to 20 programs





2900 VERTICAL AUTOCLAVE



- Three times pre-vacuum
- High quality stainless steel material
- Microcomputer automatic control, arbitrary setting sterilization parameters
- Safety interlocking device
- LED touch screen display
- Standard test interface
- Self-expanding seal
- With drying function
- With stainless steel bucket or basket
- With automatic protection function: over temperature protection; over pressure self discharge protection; low water level protection; anti dry burning.
- Automatic water
- External steam Generator
- Buzzer reminder after sterilization, automatic stop.
- Automatic discharge of cold air, automatic steam exhaust after sterilization
- Internal printer

SPECIFICATIONS

Model	BAVT-2901	BAVT-2902	BAVT-2903	BAVT-2904
Capacity	35 L	50 L	75 L	100 L
Power	3.5 kW		4.5 kW	
Voltage	AC 220V. 50HZ			
Design pressure	0.25 Mpa			
Design temperature	139°C			
Rated working pressure	0.22 Mpa			
Rated working temperature	134°C			
Sterilization temperature setting range	116°C-134°C			
Sterilization time setting range	4 ~ 120 min			
Drying time setting range	0 ~ 240 min			
Chamber volume	Ø 350x400	Ø 400x505	Ø 400x655	Ø 450x650
Bucket dimension	Ø 330x320	Ø 380x430	Ø 380x560	Ø 420x540
Basket dimension	Ø 320x350	Ø 360x280x1	Ø 360x280x2	Ø 410x300x2
Outer dimension	698x498x940	630x500x1180	725x525x1170	760x560x1260
Packing size	790x610x1130	800x600x1230	840x610x1310	870x650x1380
Net weight	115 Kg	118 Kg	124 Kg	132 Kg
Gross Weight	130 Kg	138 Kg	150 Kg	159 Kg



BAVT-2901



BAVT-2902



BAVT-2903



BAVT-2904

3000 VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure
 High quality stainless steel material
 With drying function
 Microcomputer automatic control, arbitrary setting sterilization parameters
 Safety interlocking device
 Digital display and Touch buttons
 Standard test interface
 Self-expanding seal
 With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning
 With stainless steel bucket or basket
 Buzzer reminder after sterilization, automatic stop
 Automatic discharge of cold air, automatic steam exhaust after sterilization, Internal steam cycle

SPECIFICATIONS

Model	BAVT-3001	BAVT-3002	BAVT-3003	BAVT-3004
Capacity	35 L	50 L	75 L	100 L
Power	3.5 kW		4.5 kW	
Voltage	AC 220V. 50HZ			
Design pressure	0.25 Mpa			
Design temperature	139℃			
Rated working pressure	0.22 Mpa			
Rated working temperature	134℃			
Sterilization temperature setting range	116℃-134℃			
Sterilization time setting range	4 ~ 120 min			
Drying time setting range	0 ~ 240 min			
Chamber volume	Ø 350x400	Ø 350x525	Ø 400x625	Ø 450x650
Bucket dimension	Ø 330x320	Ø 330x460	Ø 380x560	Ø 420x540
Basket dimension	Ø 320x350	Ø 320x240x2	Ø 360x280x2	Ø 410x300x2
Outer dimension	698x498x940	630x500x1180	725x525x1170	788x588x1255
Packing size	790x610x1130	800x600x1230	840x610x1310	870x650x1380
Net weight	110 Kg	112 Kg	114 Kg	129 Kg
Gross Weight	125 kg	130 Kg	145 Kg	156 Kg



BAVT-3001



BAVT-3002



BAVT-3003



BAVT-3004

3100 VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure

High quality stainless steel material

Microcomputer automatic control, arbitrary setting sterilization parameters

Yr Safety interlocking device

Digital display and Touch buttons

Standard test interface

Self-expanding seal

With drying function

With automatic protection function: over temperature protection; over pressure protection; low water level protection; anti dry burning

With stainless steel bucket or basket Buzzer reminder after sterilization, automatic stop

Automatic discharge of cold air, automatic steam exhaust after sterilization

SPECIFICATIONS

Model	BAVT-3101	BAVT-3102	BAVT-3103	BAVT-3104
Capacity	35 L	50 L	75 L	100 L
Power	3.5 kW		4.5 kW	
Voltage	AC 220V. 50HZ			
Design pressure	0.25 Mpa			
Design temperature	139°C			
Rated working pressure	0.22 Mpa			
Rated working temperature	134°C			
Sterilization temperature setting range	116°C-134°C			
Sterilization time setting range	4 ~ 120 min			
Drying time setting range	0 ~ 240 min			
Chamber volume	Ø 350x400	Ø 350x525	Ø 400x625	Ø 450x650
Bucket dimension	Ø 330x320	Ø 330x460	Ø 380x560	Ø 420x540
Basket dimension	Ø 320x350	Ø 320x240x2	Ø 360x280x2	Ø 410x300x2
Outer dimension	480x480x1000	480x480x1120	525x525x1190	550x550x1250
Packing size	580x630x1130	580x630x1230	630x680x1290	690x650x1410
Net weight	58 Kg	80 Kg	90 Kg	106 Kg
Gross Weight	79 Kg	106 Kg	120 Kg	133 Kg



BAVT-3101



BAVT-3102



BAVT-3103



BAVT-3104

3200 VERTICAL AUTOCLAVE



Hand wheel translation and quick opening structure

High quality stainless steel material

Microcomputer automatic control, arbitrary setting sterilization parameters

Safety interlocking device

Digital display and Touch buttons

Standard test interface

Self-expanding seal

With automatic protection function: over temperature protection; over pressure self discharge protection; low water level protection; anti dry burning

With stainless steel bucket or basket

Buzzer reminder after sterilization, automatic stop

Automatic discharge of cold air, automatic steam exhaust after sterilization

SPECIFICATIONS

Model	BAVT-3201	BAVT-3202	BAVT-3203	BAVT-3204
Capacity	35 L	50 L	75 L	100 L
Power	3.5 kW		4.5 kW	
Voltage	AC 220V. 50HZ			
Design pressure	0.25 Mpa			
Design temperature	139℃			
Rated working pressure	0.22 Mpa			
Rated working temperature	134℃			
Sterilization temperature setting range	116℃-134℃			
Sterilization time setting range	4 ~ 120 min			
Chamber volume	Ø 350x400	Ø 350x525	Ø 400x625	Ø 450x650
Bucket dimension	Ø 330x320	Ø 330x480	Ø 380x540	Ø 430x550
Basket dimension	Ø 320x350	Ø 320x240x2	Ø 360x280x2	Ø 410x300x2
Outer dimension	480x480x1000	480x480x1120	525x525x1190	555x555x1250
Packing size	580x630x1130	580x630x1230	630x680x1290	690x650x1410
Net weight	58 Kg	80 Kg	90 Kg	106 Kg
Gross Weight	79 Kg	106 Kg	120 Kg	133 Kg



BAVT-3201



BAVT-3202



BAVT-3203



BAVT-3204

3300 VERTICAL AUTOCLAVE



Bolt structure

High quality stainless steel material

Microcomputer automatic control, arbitrary set of sterilization parameters

Digital display and touch buttons

Standard test interface

The container with automatic protection function: over temperature protection; over pressure protection self discharge; low water level protection, anti dry burning

With stainless steel bucket or basket

Buzzer reminder after sterilization, automatic stop.

Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BAVT-3301	BAVT-3302	BAVT-3303
Capacity	35 L	50 L	75 L
Power	2.5 kW	3.5 kW	4.5 kW
Voltage	AC 220V. 50HZ		
Design pressure	0.25 Mpa		
Design temperature	139°C		
Rated working pressure	0.22 Mpa		
Rated working temperature	134°C		
Sterilization temperature setting range	116°C-134°C		
Sterilization time setting range	4 ~ 120 min		
Chamber volume	Ø 300x500	Ø 350x525	Ø 400x625
Bucket dimension	Ø 280x400	Ø 330x460	Ø 380x560
Outer dimension	400x400x860	480x480x960	525x525x1060
Packing size	550x500x1030	580x630x1100	630x680x1190
Net weight	41 Kg	53 Kg	61 Kg
Gross Weight	58 Kg	75 Kg	82 kg
Basket dimension	-	Ø 320x240x2	Ø 360x280x2

Model	BAVT-3304	BAVT-3305	BAVT-3306
Capacity	100 L	120 L	150 L
Power	4.5 kW		6.0 kW
Voltage	AC 220V. 50HZ		
Design pressure	0.25 Mpa		

Design temperature	139°C		
Rated working pressure	0.22 Mpa		
Rated working temperature	134°C		
Sterilization temperature setting range	116°C-134°C		
Sterilization time setting range	4 ~ 120 min		
Chamber volume	Ø 450x650	Ø 450x760	Ø 500x780
Bucket dimension	Ø 420x540	Ø 420x660	Ø 470x720
Outer dimension	555x555x1080	555x599x1200	605x649x1210
Packing size	690x650x1230	690x650x1330	750x700x1370
Net weight	78 Kg	82 Kg	90 Kg
Gross Weight	103 Kg	108 Kg	116 Kg
Basket dimension	Ø 410x300x2	Ø 410x350x2	Ø 460x350x2



BAVT-3301



BAVT-3302



BAVT-3303



BAVT-3304



BAVT-3305



BAVT-3306



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com