



AIR SHOWER PASS BOX

AIR SHOWER PASS BOX

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

100 AIR SHOWER PASS BOX



Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BAPB-101	BAPB-102	BAPB-103
Clean Air Speed	0.36~0.54m/s		
UV Lamp	15Wx1		
External Size (WxDxH)	660x570x1000 mm	760x670x1100 mm	860x770x1200 mm
Internal Size (WxDxH)	500x500x500 mm	600x600x600 mm	700x700x700 mm
Material	304 stainless steel		
Interlock	Electronical		
HEPA Filter	Efficiency ≥99.995% at 0.3μm		
Gross Weight	80 kg	100 kg	120 kg
Packing Size(WxDxH)	780x700x1250 mm	880x800x1350mm	980x900x1450mm
Power Supply	220V, 50Hz; 110V, 50/60Hz		



BAPB-101



BAPB-102



BAPB-103



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com