



PASS BOX

## PASS BOX

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.

### 100 PASS BOX



Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

## SPECIFICATIONS

| Model                 | BSPB-101                  | BSPB-102       | BSPB-103       |
|-----------------------|---------------------------|----------------|----------------|
| UV Lamp               | 8Wx1                      |                |                |
| External Size (WxDxH) | 660x570x580 mm            | 760x670x680 mm | 860x770x780 mm |
| Internal Size (WxDxH) | 500x500x500 mm            | 600x600x600 mm | 700x700x700 mm |
| Material              | 304 stainless steel       |                |                |
| Interlock             | Electronical              |                |                |
| Gross Weight          | 55 kg                     | 70kg           | 90kg           |
| Packing Size(WxDxH)   | 780x700x770 mm            | 880x800x870 mm | 980x900x970 mm |
| Power Supply          | 220V, 50Hz; 110V, 50/60Hz |                |                |



BSPB-101



BSPB-102



BSPB-103



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)