



AIR SHOWER

AIR SHOWER

Used in Laboratory, Microbiological, Research ,University ,,
Also known as Personal Cleaning Booth.

100 AIR SHOWER



Photoelectric sensor, automatic blowing

Adopt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time (0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BASB-101	BASB-102	BASB-103
External Size (WxDxH)	1300x1000x2050 mm	1500x1000x2050 mm	1500x2000x2050 mm
Internal Size (WxDxH)	800x900x1930 mm		800x1900x1930 mm
HEPA filter	99.995 % efficiency at 0.3 μ m		
Spraying Wind Speed	25 m/s		
Air Shower time	0-99 s (adjustable)		
No. of SS Nozzels	6	12	24
Fluorescent Lamp	12Wx1		
Consumption	550 W	1100 W	2200 W
Cabinet Material (Optional)	1.Full cold-rolled steel with powder coating 2.Internal stainless steel, external with cold-rolled steel (optional) 3.Full stainless steel (optional)		
Door Material (Optional)	1.Cold-rolled steel with powder coating 2.Stainless steel (optional)		
Chassis	Stainless Steel		
Power Supply	AC220 V \pm 10%,50/60 Hz ; 110 V \pm 10%, 60 Hz		
Gross Weight	300 kg	400 kg	700 kg
Package Size (WxDxH)	1400x1200x2200 mm	1600x1200x2200 mm	1600x2200x2200 mm



BASB-101



BASB-102



BASB-103



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com