



LABORATORY HIGH SPEED CENTRIFUGE

LABORATORY HIGH SPEED CENTRIFUGE

Designed around your applications to provide you competent sample processing and reliable results. Microcomputer programmable with excellent temperature controls and low noise operation maximize your productivity. Highly efficient with low maintenance requirements, it is an ideal separation tool for multiple research applications.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

BCBH-102 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC motor, no carbon dust pollution, unnecessary to maintain

Microprocessor control, small LCD display in the speed, time, RCF in operation, speed raising and reducing quick, operate simply. Running parameters can be edited

Automatically electric lid lock, super speed, imbalance protection. The centrifuge body is made of high quality steel, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction

It is with CE & ISO9001 & ISO13485 Certificates

SPECIFICATIONS

Model	BCBH-102
Maximum Capacity (No of tubes x Vol.)	10x5 ml
Maximum Speed	16000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	19040xg
Time Range	0-99 min
Overall Dimension	335x270x190 mm
Noise Level	<65 dB(A)
Net. Weight	15 kg
Power Supply	AC220V 50Hz/110V 60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2300807006	Angle rotor	40x0.2 ml	16000	19040xg
2300807007	Angle rotor	24x0.5 ml	16000	18480xg
2300807008	Angle rotor	12x1.5/2 ml	16000	16260xg
2300807009	Angle rotor	10x5 ml	16000	17880xg
2300807010	Angle rotor	20x1.5/2 ml	14000	15580xg

200 BENCHTOP HIGH SPEED CENTRIFUGE



Microcomputer controlled programmable operation with touch panel

TFT true-colour LCD wide screen displays set and run conditions for easy monitoring

Brushless DC Motor ensures low maintenance, high efficiency, constant speed of rotation

Store 20 programs in centrifuge memory which are retained even after powering down the centrifuge

Fully stainless steel structure and a stainless steel cavity makes the machine safe and more efficient

Acceleration/Deceleration takes around 30 sec to 1 min

Multiple layers of shock absorbers with auto balancing function

Dedicated key indicating centrifugal force

System is complemented by multiple new accessories, including tube racks, centrifuge adapters and rotors

Safety features: Self diagnostic, Rotor unbalance detector, Lid protection, Over speed protection, Over temperature protection

SPECIFICATIONS

Model	BCBH-201	BCBH-202	BCBH-204
Maximum Speed	12000 rpm	16500 rpm	18500 rpm
Speed Increment	±50 rpm		±30 rpm
Maximum RCF	13800xg	19600xg	24640xg
Time Range	0-99h 59 min	1min-99h59min59s	99h59min59s,
Overall Dimension	500x440x330 mm	500mm×470mm×330mm	500x470x330 mm
Program	20 user programs		
Control	Microcomputer		
Display	LCD		
Noise Level	<65dB(A)		
Weight (Net/Gross)	35/42 kg	18/30 kg	35/42 kg
Power	370 W	350 W	500 W
Power Supply	220V, 50Hz	AC220V50Hz	
Maximum Capacity (No of tubes x Vol.)	-	4x100 ml	4x100ml

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}	For Models
2300909006	Fixed Angle Rotor	12x1.5/2.0 ml	18500	23610xg	BCBH-204
2300909007	Fixed Angle Rotor	18x1.5/0.5 ml	15000	17770xg	BCBH-204
2300909008	Fixed Angle Rotor	24x1.5/2.0 ml	14000	18180xg	BCBH-204
2300909009	Fixed Angle Rotor	10x5 ml	14000	13520xg	BCBH-204
2300909010	Fixed Angle Rotor	8x15 ml	12000	14800xg	BCBH-204
2300909011	Fixed Angle Rotor	12x10ml	12000	13400xg	BCBH-204
2300909012	Fixed Angle Rotor	12x5ml	15000	15600xg	BCBH-204
2300909013	Fixed Angle Rotor	6x50 ml	11000	11610xg	BCBH-204
2300909014	Fixed Angle Rotor	6x50 ml	6000	3460xg	BCBH-204
2300909015	Swinging-bucket Rotor	12x10 ml	4000	2330xg	BCBH-204

2300909016	Swing Micro Plate Rotor	2x2x48 holes	3000	1400xg	BCBH-204
------------	-------------------------	--------------	------	--------	----------



BCBH-201



BCBH-202



BCBH-204

BCBH-203 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running

Automatically electric lid lock, super speed, over temperature protection and imbalance protection

The centrifuge body is made of high quality steel, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable

Microprocessor control, TD5A with digital display and TD5C with LCD display which indicate the speed, time, RCF in operation, speed raising and reducing quick, operate simply

10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice

Three layers protection of steel structure and get the ideal centrifugation result

SPECIFICATIONS


Model	BCBH-203
Maximum Capacity (No of tubes x Vol.)	6x100 ml
Maximum Speed	16000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	20600xg
Time Range	1~99 min
Overall Dimension	513x370x320 mm
Noise Level	<60 dB(A)
Net. Weight	42 kg
Power Supply	AC 220V 50HZ 10A
Package	carton box

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2300908006	Angle rotor	12x1.5/2 ml	16000	16270xg
2300908007	Angle rotor	40x0.5 ml	14000	19970xg
2300908008	Angle rotor	24x1.5/2 ml	15000	20600xg
2300908009	Angle rotor	30x1.5/2 ml	13500	19340xg
2300908010	Angle rotor	16x5 ml	15000	19350xg

2300908011	Angle rotor	12x7 ml	14000	16370xg
2300908012	Angle rotor	12x15 ml	10000	11840xg
2300908013	Angle rotor	12x10 ml	12000	14510xg
2300908014	Angle rotor	8x20 ml	12000	14510xg
2300908015	Angle rotor	6x30 ml	12000	14000xg
2300908016	Angle rotor	6x50 ml	10000	11140xg
2300908017	Angle rotor	6x70 ml	10000	10810xg
2300908018	Angle rotor	4x100 ml	12000	14850xg
2300908019	Angle rotor	6x100 ml	8000	7280xg
2300908020	Angle rotor	6x10 ml	14000	16460xg
2300908021	Angle rotor	30x0.5 ml	15000	18510xg
2300908022	Angle rotor	24pieces capillary vessel	12000	15800xg
2300908023	Vertical rotor	16x5ml ml	16000	16540xg
2300908024	Bucket rotor	2x3x48 well	4000	2300xg

BCBH-302 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running

Automatically electric lid lock, super speed, over temperature protection and imbalance protection

The centrifuge body is made of high quality steel, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable

Microprocessor control, TD5A with digital display and TD5C with LCD display which indicate the speed, time, RCF in operation, speed raising and reducing quick, operate simply

10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice

Three layers protection of steel structure and get the ideal centrifugation result

SPECIFICATIONS

Model	BCBH-302
Maximum Capacity (No of tubes x Vol.)	6x100 ml
Maximum Speed	21000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	28030xg
Time Range	1~99 min
Overall Dimension	513x370x320 mm
Noise Level	<60 dB(A)
Net. Weight	42 kg
Power Supply	AC 220V 50HZ 10A
Package	carton box

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301008006	Angle rotor	12x1.5/2 ml	21000	28030 _{xg}
2301008007	Angle rotor	24x1.5/2 ml	17000	26460 _{xg}
2301008008	Angle rotor	16x5 ml	16000	22020 _{xg}
2301008009	Angle rotor	6x30 ml	14000	19060 _{xg}
2301008010	Angle rotor	12x10 ml	15000	22680 _{xg}
2301008011	Angle rotor	12x15 ml	10000	11840 _{xg}
2301008012	Angle rotor	6x10 ml	16000	21500 _{xg}
2301008013	Angle rotor	6x50 ml	12000	16050 _{xg}
2301008014	Angle rotor	8x20 ml	15000	22680 _{xg}
2301008015	Angle rotor	6x70 ml	12000	15570 _{xg}
2301008016	Angle rotor	6x100 ml	10000	11380 _{xg}
2300908024	Bucket rotor	2x3x48 well	4000	2300 _{xg}
2301008017	Angle rotor	40x0.5 ml	15000	22920 _{xg}
2301008018	Angle rotor	30x1.5/2 ml	13500	19340 _{xg}
2301008019	Angle rotor	12x7 ml	16000	21380 _{xg}
2301008020	Angle rotor	4x100 ml	12000	14850 _{xg}
2301008021	Angle rotor	30x0.5 ml	18000	26660 _{xg}
2301008022	Angle rotor	48x1.5/2 ml	13000	17930 _{xg}
2301008023	Vertical rotor	16x5 ml	14000	12660 _{xg}

BCBH-304 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless motor, no pollution, maintenance-free

Microprocessor control, LCD display which indicates the speed, time, RCF in operation, 10 kinds of brake setting, operate simply

Electric lid lock, super speed, imbalance protection

Centrifuge body is made of high quality steel, stainless steel chamber, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction

3 tiers protection steel cover and get the ideal centrifugation result

SPECIFICATIONS

Model	BCBH-304
Maximum Capacity (No of tubes x Vol.)	6x10 ml
Maximum Speed	16000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	19040 _{xg}
Time Range	0-99 min
RPM / RCF convert	Yes
Acceleration/ Deceleration	10 kinds

Overall Dimension	370x290x215 mm
Control	Microcomputer
Noise Level	<65 dB(A)
Net. Weight	16 kg
Power	470 W
Power Supply	AC 200 50 Hz / 110 V 60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301010006	Angle rotor	40x0.2 ml	16000	19040xg
2301010007	Angle rotor	12x1.5/2 ml	16000	16260xg
2301010008	Angle rotor	10x5 ml	16000	17880xg
2301010009	Angle rotor	24x1.5/2 ml	14000	17950xg
2301010010	Angle rotor	24x0.5 ml	16000	18480xg
2301010011	Angle rotor	4×8PCR	14000	12070xg
2301010012	Angle rotor	6×10ml	13000	14190xg
2301010013	Angle rotor	24 pieces capillary vessel	12000	15800xg
2301010014	Angle rotor	40x0.5ml	13000	17210xg
2301010015	Angle rotor	30x0.5ml	13000	13900xg

BCBH-305 BENCHTOP HIGH SPEED CENTRIFUGE



10 kinds of acceleration and deceleration

Multiple rotor for your choice

High speed, low noise, compact design and small vibration

User-friendly LED display provides speed, time and temperature parameters

Brushless frequency motor which in great torque, no powder pollution, free maintenance, quick in speed up and down

Imbalance sensor and over speed sensor, the centrifuge body is made of high-quality steel, very safe

SPECIFICATIONS

Model	BCBH-305
Maximum Capacity (No of tubes x Vol.)	4x300ml
Maximum Speed	16000 rpm
Speed Accuracy	±20 rpm
Maximum RCF	27170xg
Time Range	0-99h 59 min
RPM / RCF convert	Yes
Acceleration/ Deceleration	10 kinds
Control	Microcomputer
Noise Level	<55 dB(A)

Power Supply	AC 200 50 Hz / 110 V 60Hz
--------------	---------------------------

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301011006	Swing rotor rack (aluminium alloy)			
2301011007	Swing buckets (round cup) (aluminium alloy)	4x300 ml	4500	4030xg
2301011008	Swing buckets (square cup) (aluminium alloy)	4x300 ml	4500	3990xg
2301011009	Swing buckets (aluminium alloy)	4x24x7 ml	4000	3150xg
2301011010	Swing buckets (aluminium alloy)	4x24x5 ml	4000	2750xg
2301011011	Swing buckets (aluminium alloy)	4x16x10 ml	4000	3150xg
2301011012	Microplate rotor	2x4x96 well	4000	2910xg
2301011013	Swing rotor rack (stainless steel)		4000	2910xg
2301011014	Swing buckets (stainless steel)	4x1x50 ml	5000	4730xg
2301011015	Swing buckets (stainless steel)	4x1x100 ml	5000	4730xg
2301011016	Swing buckets (stainless steel)	4x2x100 ml	4000	3020xg
2301011017	Swing buckets (stainless steel)	4x2x50 ml	4000	3020xg
2301011018	Swing buckets (stainless steel)	4x8x15 ml	4000	3200xg
2301011019	Swing buckets (stainless steel)	4x6x15 ml	4000	3200xg
2301011020	Swing buckets (stainless steel)	4x4x10 ml	4000	2770xg
2301011021	Swing buckets (stainless steel)	4x4x7/5 ml vacuum tube adapter	4000	2360xg
2301011022	Swing buckets (stainless steel)	4x6x10 ml	4000	2770xg
2301011023	Swing buckets (stainless steel)	4x6x7/5 ml vacuum tube adapter	4000	2360xg
2301011024	Swing buckets (stainless steel)	4x8x10 ml	4000	2770xg
2301011025	Swing buckets (stainless steel)	4x8x7/5 ml vacuum tube adapter	4000	2360xg
2301011026	Swing buckets (stainless steel)	4x16x5 ml vacuum tube	4000	2520xg
2301011027	Swing buckets (stainless steel)	4x16x7 ml vacuum tube	4000	2880xg
2301011028	Swing buckets (stainless steel)	4x12x10 ml vacuum tube	4000	2880xg
2301011029	Fixed Rotor	4x8PCR	16000	15760xg
2301011030	Fixed Rotor	6x8PCR	16000	24370xg
2301011031	Fixed Rotor	8x8PCR	16000	17480xg
2301011032	Fixed Rotor	12x8PCR	16000	26090xg
2301011033	Fixed Rotor	40x0.5 ml	15000	22920xg
2301011034	Fixed Rotor	12x1.5/2 ml	16000	17940xg
2301011035	Fixed Rotor	24x1.5/2 ml	16000	23440xg
2301011036	Fixed Rotor	30x1.5/2 ml	16000	27170xg
2301011037	Fixed Rotor	48x1.5/2 ml	16000	27170xg
2301011038	Fixed Rotor	16x5 ml	16000	22020xg
2301011039	Fixed Rotor	6x10 ml	16000	21500xg
2301011040	Fixed Rotor	12x10 ml	16000	25800xg
2301011041	Fixed Rotor	16x10 ml	14000	22610xg
2301011042	Fixed Rotor	8x15 ml	14000	20630xg
2301011043	Fixed Rotor	12x15 ml	12000	17060xg
2301011044	Fixed Rotor	24x15 ml	5000	3500xg
2301011045	Fixed Rotor	30x15 ml	5000	3830xg
2301011046	Fixed Rotor	6x30 ml	15000	21880xg
2301011047	Fixed Rotor	12x30 ml	12000	17640xg

2301011048	Fixed Rotor	6x50 ml (falcon)	14000	21850xg
2301011049	Fixed Rotor	8x50 ml (falcon)	12000	18200xg
2301011050	Fixed Rotor	10x50 ml (falcon)	9000	11240xg
2301011051	Fixed Rotor	6x50 ml (round)	14000	21730xg
2301011052	Fixed Rotor	8x50 ml (round)	12000	17880xg
2301011053	Fixed Rotor	4x85 ml	14000	21960xg
2301011054	Fixed Rotor	6x85 ml	12000	17750xg
2301011055	Fixed Rotor	8x85 ml	9000	11160xg
2301011056	Fixed Rotor	4x100 ml	14000	20220xg
2301011057	Fixed Rotor	6x100 ml	12000	16390xg
2301011058	Fixed Rotor	8x100 ml	9000	11120xg
2301011059	Fixed Rotor	12x100 ml	4000	2970xg
2301011060	Fixed Rotor	4x250 ml	8000	9090xg
2301011061	Fixed Rotor	24 capillaries	14000	15800xg
2301011062	Swing Rotor	4x5 ml	15000	19920xg
2301011063	Swing Rotor	4x30 ml	12000	20090xg
2301011064	Swing Rotor	4x50 ml	10000	16240xg
2301011065	Vertical rotor	16x5 ml	16000	16450xg
2301011066	Vertical rotor	8x30 ml	14000	19750xg

BCBH-306 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC frequency motor with simpler construction, more reliable performance, longer life and quietly running

Microprocessor control, LED display in speed, temperature and time etc.

Automatically electric lid lock, super speed, over temperature protection and imbalance protection

The centrifuge body is made of high quality steel, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction, safe and reliable

There are several kinds of rotors for your choice, adapters can be designed by experimental requirements

3 layers protection of steel structure, safe and reliable

SPECIFICATIONS

Model	BCBH-306
Maximum Capacity (No of tubes x Vol.)	4x800ml
Maximum Speed	21000
Speed Accuracy	±20 rpm
Maximum RCF	30910xg
Time Range	0-99h 59 min
Overall Dimension	700x500x400 mm
Control	Microcomputer
Noise Level	<58 dB(A)
Net. Weight	75 kg
Power Supply	AC110/ 220V 50/60HZ

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	RPM	RCF _{xg}
2301012006	Angle rotor	12x1.5/2ml	21000	28030xg
2301012007	Angle rotor	40x0.5ml	15000	22920xg
2301012008	Angle rotor	24x1.5/2ml	17000	26460xg
2301012009	Angle rotor	30x1.5/2ml	13500	19340xg
2301012010	Angle rotor	16x5ml	16000	22020xg
2301012011	Angle rotor	12x7ml	16000	21380xg
2301012012	Angle rotor	12x15ml	10000	11840xg
2301012013	Angle rotor	12x10ml	15000	22680xg
2301012014	Angle rotor	8x20ml	15000	22680xg
2301012015	Angle rotor	6x30ml	14000	19060xg
2301012016	Angle rotor	6x50ml	12000	16050xg
2301012017	Angle rotor	6x70ml	12000	15570xg
2301012018	Angle rotor	4x100ml	12000	14850xg
2301012019	Angle rotor	6x100ml	10000	11380xg
2301012020	Angle rotor	6x10ml	16000	21500xg
2301012021	Angle rotor	30x0.5ml	18000	26660xg
2301012022	Angle rotor	48x1.5/2ml	13000	17930xg
2301012023	Bucket rotor	4x4x96 well	4000	2840xg
2301012024	Swing rotor	4x5 ml	15000	19920xg
2301012025	Bucket rotor	2x4x96 well	4000	2490xg
2301012026	Swing rotor	4x40x5 ml	4000	2890xg
2301012027	Swing rotor	4x40x7 ml	4000	3270xg
2301012028	Swing rotor	4x28x10 ml	4000	3130xg
2301012029	Vertical rotor	16x5 ml	16000	16450xg
2301012030	Swing rotor	4x800 ml	4000	3450xg
2300908024	Bucket rotor	2x3x48 well	4000	2300xg

BCBH-307 BENCHTOP HIGH SPEED CENTRIFUGE



Brushless DC motor, no carbon dust pollution, free maintenance

Microprocessor control, small LCD display

Automatic rotor identification function

Electronic lock, over speed and imbalance protection

10 kinds of program stored in the memory, 10 kinds of accelerating and decelerating speed for your choice

It is with CE & ISO9001 & ISO13485 Certificates

SPECIFICATIONS

Model	BCBH-307
Maximum Capacity (No of tubes x Vol.)	24 capillary tubes
Maximum Speed	12000 rpm
Speed Accuracy	±20rpm
Maximum RCF	15800xg
Time Range	0~99min
Noise Level	≤55 dBA
Power Supply	AC 220 50Hz 2A
Overall Dimension	335x270x190 mm
Net. Weight	13 kg
Package	Carton box

BCBHR-211 BENCHTOP HIGH SPEED REFRIGERATED CENTRIFUGE



In static pre-cooling mode, when the centrifuge cover is closed, if the set temperature is lower than the ambient temperature, the compressor will automatically starts to cool

The sample temperature can still maintain at 4°C when running at the highest speed, intermittent temperature control mode, fluorine-free refrigeration, low carbon and environmental protection

Strong cooling capacity, the centrifuge chamber can be lowered from room temperature to 4°C in less than 8 minutes

High-speed and stable operation, the rotor and the motor shaft are matched with high-precision cones, with high concentricity and low noise

High-strength metal protective rings are used in the upper and lower shell, which have high structural safety and ensure the safe operation of the rotor at high speed

Patent of the upper cover safety lock technology ensures the safety during operation, and has one-button cover opening function

Powerful DC brushless motor, maintenance-free, fast drive, ultra-quiet operation, low noise

Soft-brake function to protect special sensitive samples, and the lifting speed is adjustable to prevent the sample from being resuspended

Using 4.3-inch color LCD screen, can display speed, centrifugal force and time at the same time

Has the function of beeping sound when operation completed, and has the function of running status display and error code prompting

There is a PROG key to call the built-in 10 groups of programs with one key, which is convenient and quick to use

Equipped with an instant centrifugation button, it can be centrifuged quickly to meet more experimental needs

High-strength aluminum alloy rotor, no service life limit, can be sterilized at high temperature and high pressure

One machine multi-purpose, suitable for 0.2, 0.5, 1.5, 2.0, 5.0ml centrifuge tubes and PCR tubes, and 4 rotors are optional

After the rotor is replaced, the rotor model can be automatically identified, and no need to manually set the rotor model

SPECIFICATIONS

Model	BCBHR-211
Rotor capacity	0.2 ml/0.5 ml/1.5 ml/2 ml x 24; 0.2 ml x 8 strip x 4; 0.5 ml x 36; 5 ml x 18;
Speed Range	500 - 15000
Maximum RCF	21400 x g
Speed up down time range	15 or 45 s (to protect special sensitive samples from resuspension)
Time range	1 s~99 min 59 s (can be centrifuged instantly)
Other functions	Speed/acceleration display function, instantaneous centrifugal function, running progress display, sound prompt function
Weight (Net/Gross)	26 Kg
Overall Dimension	W.335 X D.574 X H.234 mm 26 kg
Voltage	AC 220 V 50 Hz
Working noise	≤ 57 dB
Safety performance	High-strength metal protective ring, safety cover lock, overspeed, overtemperature, status diagnosis
Temp. range	- 20°C to + 40°C
Temp.control accuracy (at 4°C)	± 2 °C
Time from R. T. to 4°C	≤ 8 min
Fuse	250 V ,15 A, φ 5 × 20
Power Supply	1000 W

BCBH-401 BENCHTOP HIGH SPEED CENTRIFUGE

Widely used in the field of biochemistry, biological products, pharmaceutical factory and laboratory

Brushless frequency motor for model G12 which in great torque, free maintenance, no powder pollution, quick in speed up and down

Model G12 with the digital display which indicates the speed、time

Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice

There are high speed rotors and low speed rotors for your choice

Electric lid lock, compact design, super speed and imbalance protection

7The centrifuge body is made of high-quality steel, safe and reliable

SPECIFICATIONS

Model	BCBH-401
Maximum Capacity (No of tubes x Vol.)	4 x 1000 ml
Maximum Speed	12000 rpm
Speed Accuracy	± 20 r / min
Maximum RCF	22650 × g
Time Range	1 ~ 9 h 59 min
Acc / Dec	10 Kinds
Noise Level	≤ 60 dB(A)
Certificates	CE, ISO & Calibration report are available
RPM / RCF Convert	Yes
Overall Dimension	780 × 675 × 830 mm

Net. Weight	200 kg
Power Supply	AC 220V / 110 V 50 HZ / 60 HZ

BCBH-402 BENCHTOP HIGH SPEED CENTRIFUGE

Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down

Microprocessor control

Digital display in time, speed and RCF

10 kinds of accelerating and decelerating speed for your choice

Automatically electric lid lock, over speed, over temperature protection and imbalance protection

High quality steel centrifuge body, stainless steel chamber, 3 tiers protection, safe and reliable

Rotor is connected to spindle by specialized taper sleeve, loading simple and quick, no direction

SPECIFICATIONS

Model	BCBH-402
Maximum Capacity (No of tubes x Vol.)	6 x 100 ml
Maximum Speed	18500 rpm
Speed Accuracy	± 20 r / min
Maximum RCF	23990 × g
Time Range	0 ~ 99 min
Acc / Dec	10 Kinds
Noise Level	≤ 65 dB(A)
Certificates	CE, ISO & Calibration report are available
RPM / RCF Convert	Yes
Overall Dimension	513 × 370 × 320 mm
Net. Weight	42 kg
Power Supply	AC 220 V / 110 V 50 HZ / 60 HZ



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com