



BULLET ICE MAKER

BULLET ICE MAKER

Biolab Bullet Ice Maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of bullet shaped ice cubes with ice making capacity from 20 kg/24h up to 120 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

100 BULLET ICE MAKER



- Vertical stirring engine model
- Refrigeration system produces uniform shape ice bullets
- Bell-shaped ice with a hexagon hole
- Highly durable stainless steel 304 2B built-up
- Automatic controller increases working reliability
- Large ice storage bin
- Eco-friendly, energy saving non-Freon refrigerant

SPECIFICATIONS

Model	BIBU-101	BIBU-102	BIBU-104	BIBU-109
Capacity	20 kg/24 h	25 kg/24 h	50 kg/24 h	80 kg/24 h
Storage Bin Capacity	6 kg	7 kg	14 kg	25 kg
Refrigerant	R134a			
Condenser Type	Air Cooled			
Overall Dimension	330x470x605 mm	400x543x670 mm	400x510x825 mm	500x612x900 mm
Packing Dimension	405x555x665 mm	470x615x730 mm	470x590x875 mm	560x670x950 mm
Weight	28/32 kg	30/33 kg	35/38 kg	51/56 kg
Power	220 W	240 W	380 W	580 W
Power Supply	220V /110 V, 50 / 60 Hz			

OPTIONAL ACCESSORIES

Accessory Code	Name
2700606006	Water filter



BIBU-101

BIBU-102

BIBU-104

BIBU-109



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com