



ICE FLAKER

ICE FLAKER

Biolab Flake Ice maker is compact in size, made of highly durable material with energy efficient features. Increased efficiency gives quick production of ice flakes with ice making capacity from 50 kg/24h up to 130 kg/24h. Self contained storage bin eliminates the need of a separate storage unit.

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

200 ICE FLAKER



- High quality stainless steel construction
- Fully automatic operation system
- Fluorine-free foam thermal cover ensures good insulation
- High speed compact clear flake ice-production
- Advanced self-checking mode function
- Power-overloading, water shortage and ice storage full protection
- Large output and energy efficient
- Easy to clean storage bin

SPECIFICATIONS

Model	BIFL-201	BIFL-202	BIFL-203	BIFL-204
Capacity	20 kg/24 h	30 kg/24 h	40 kg/24 h	50 kg/24 h
Storage Bin Capacity	10 kg		15 kg	
Refrigerant	R134a			
Condenser Type	Air cooled			
Overall Dimension	330x470x605 mm		400x543x720 mm	
Packing Dimension	405x555x665 mm		470x615x780 mm	
Weight	28/31 kg		36/39 kg	
Power	280 W	300 W	360 W	380 W
Power Supply	220 V, 60 Hz			

Model	BIFL-205	BIFL-206	BIFL-207	BIFL-208
Capacity	60 kg/24 h	70 kg/24 h	85 kg/24 h	100 kg/24 h
Storage Bin Capacity	25 kg			40 kg
Refrigerant	R134a			
Condenser Type	Air cooled			
Overall Dimension	400x510x825 mm			500x611x940 mm
Packing Dimension	470x590x875 mm			570x685x1000 mm
Weight	42/46 kg		45/49 kg	58/63 kg
Power	420 W	460 W	480 W	520 W
Power Supply	220 V, 60 Hz			

Model	BIFL-209	BIFL-210	BIFL-211
Capacity	120 kg/24 h	200 kg/24 h	300 kg/24 h

Storage Bin Capacity	40 kg	60 kg	150 kg
Refrigerant	R134a		
Condenser Type	Air cooled		
Overall Dimension	500x611x940 mm	680x611x1100 mm	560x800x1700 mm
Packing Dimension	570x685x1000 mm	750x675x1150 mm	650x850x1750 mm
Weight	58/63 kg	80/85 kg	95/105 kg
Power	550 W	1040 W	1360 W
Power Supply	220 V, 60 Hz		

OPTIONAL ACCESSORIES

Accessory Code	Name	For Models
2700606006	Water filter	BIFL-201, BIFL-202, BIFL-203, BIFL-204, BIFL-205, BIFL-206, BIFL-207, BIFL-208



BIFL-201



BIFL-202



BIFL-203



BIFL-204



BIFL-205



BIFL-206



BIFL-207



BIFL-208



BIFL-209



BIFL-210



BIFL-211



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com