



BIOLOGICAL INDICATOR INCUBATOR

BIOLOGICAL INDICATOR INCUBATOR

Biological indicator incubator- Series helps assure the sterility of instruments, vials, media and other products in very short time. It also provides easy evaluation of incubation temperature.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

BIBI-201 BIOLOGICAL INDICATOR INCUBATOR



Unique designed. Small size, light weight, portable, and suits on various occasions.

LCD display. Temperature and time display simultaneously.

Auto fault detect and buzzer alarm. Temperature calibration function.

24V DC power supply input, low energy consumption. Built-in over- temperature protection device, safe and reliable.

Timer function. Time range from 0 to 999 min.

Variety of standard sample blocks for choice. Convenient for replacement. Block can be customized.

SPECIFICATIONS

Model	BIBI-201
Temperature Control Range	RT+5°C~100°C
Temperature Setting Range	5°C~100°C
Temperature Stability@56°C	±0.5°C
Temp. Display Accuracy	0.1°C
Heating Speed	<6min (20°C to 56°C)
Time Range	1sec~999sec or 1min~999min
Voltage	DC 24V adapter
Power	60W
Dimension (WxDxH)	110x156x92 mm
Net Weight	1kgs

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
2905615006	Block A	15 x Φ10mm
2905615007	Block B	12 x Φ12mm
2905615008	Block C	6 x Φ15mm
2905615009	Block D	6 x Φ16mm
2905615010	Block E	18 x Φ9mm
2905615011	Customized	Customized



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com