



AIR-JACKETED MULTI GAS INCUBATOR

AIR-JACKETED MULTI GAS INCUBATOR

Used in Biomedical, Pharmaceutical, Medical, Clinical research, Laboratory.

Also known as Laboratory Air-jacketed Multi gas Incubator, Multi gas Incubator Air-jacketed.

100 AIR-JACKETED MULTI GAS INCUBATOR



151 liters internal volume
Zirconium oxide Sensor
IR CO₂ sensor
Water Level Alarm
Integrated electric siphon pump
Auto-zero
Temperature deviation (at 37°C)
Moist-heat disinfection at 90°C
Relative humidity ≥95%

SPECIFICATIONS

Model	BIMA-101	BIMA-102
Capacity	151 L	
Temperature Sensors	PT1000	
CO ₂ Control Method	Microprocessor PID	
Shelve Dimension	423x445 mm	
Internal Dimension	470Wx530Dx607H mm	
Exterior Dimension	637Wx768Dx869H mm	
Shelves	3 standard, 10 maximum	
Weight	80 kg	
Max. Load per Shelf	10 kg	
Power	≤650 VA+10%	
Power Supply	220 V±10%/ 50 Hz	



BIMA-101



BIMA-102



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com