



DELIVERY TABLE

DELIVERY TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation.

Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

BHBD-203 MANUAL OBSTETRIC BED



The main bed and auxiliary bed are connected with each other, four foot wheels are equipped on auxiliary bed and brake device is spared

Foot step hydraulic elevating, Imported Y type sealing ring and is durable

The mattress adopts high density memory sponge and is once formed, easily cleaned and sterilized

Rotary brake device which is easy for moving operating table

Base seat and upright column covering are made of high quality 304 stainless steel

SPECIFICATIONS

Model	BHBD-203
Length	1950+50 mm
Width	600+50 mm
Height Adjustment	(750~950)+50 mm
Working area length/width	650x600 mm
Backboard upward	≥50°
Backboard downward	≥10°
Seat board upwards	≥30°
Seat board downwards	≥10°
Standard accessory	A pair of arm rods, a pair of leg support, a pair of arm rods, a pair foot pedal, a stainless steel sewage basin

BHBD-204 ELECTRIC GYNECOLOGICAL TABLE



Various actions are controlled by Imported brand motor

Hand-held controller and foot controller, mutually independent and is easily for operation

Leg support frame can fold up and down electric

Unique profile design which is beautiful and elegant

Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with hidden sewage basin

The table adopts high density memory sponge and is once foam molding, with no seam

Features with one-key reset function

Super mute foot wheel and equipped with brake device

SPECIFICATIONS

Model	BHBD-204
Back board length	900 mm
Seat board length	400 mm
Table width	600 mm
Trendelenburg/rev-trend	$\geq 10^\circ/\geq 12^\circ$
Backboard Adjustment	$\geq 45^\circ$
Height adjustment (from back top to ground)	1500-1750 mm
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin, auxiliary table (already installed), a foot controller

BHBD-205 ELECTRIC OBSTETRIC BED



Imported LINAK motor various actions are controlled by remote controller

Super mute big foot wheel which is easily moved and equipped with center control breaking system.

Equipped with stainless steel sewage basin

The mattress adopts high density sponge and manufactured by high quality artificial leather without seam

With various complete accessories

Standard Accessories: A pair of leg support frame, a stainless steel sewage basin, a stand for transfusion, a remote control, a power wire

SPECIFICATIONS

Model	BHBD-205
Length	1950+50 mm
Width	800+50 mm
Height Adjustment	(550~800)+50 mm
Trendelenburg/rev-trend	$\geq 5^\circ/\geq 12^\circ$
Backboard upward	$\geq 60^\circ$

Foot board upward	≥40°
Foot board downward	≥90°
Standard accessory	A pair of leg support, a sewage basin, a power cord, a IV pole, a remote controller
Input power	150+10%
Power Supply	220+22V, 50+1Hz

BHBD-206 MANUAL MULTI-FUNCTION GYNEACOLOGICAL TABLE



Material: Steel Powder Coated S.S, Color: Blue

SPECIFICATIONS

Model	BHBD-206
Function	Multi
Back-rest lifting angle	≥20°
Back section up adjusted	≥60°
Seat board upwards	≥8°
Seat board downwards	≥10°
Load Bearing	250 kg
Height adjustment	≤650 mm
Tabletop Dimension	1240Lx600W mm
Upright Dimension	500Lx600W mm
Dimension	1750Lx600Wx650~900 mm

BHBD-208 DELIVERY TABLE 4002



Imported LINAK motor various actions are controlled by foot controller

The base seat and upright column covering are made of high quality 304 stainless steel which is anti corrosive and can be easily clean

External folding auxiliary table is for easy dismantling

Equipped with stainless steel sewage basin and other accessories

Foot treadle brake device which is easy for moving obstetric table

SPECIFICATIONS

Model	BHBD-208
Table Length	2000+50 mm
Table width	600+10 mm
Height Adjustment	(650-900)+50 mm
Trendelenburg/rev-trend	$\geq 10^\circ / \geq 12^\circ$
Backboard upward	$\geq 50^\circ$
Leg board outward	$\geq 90^\circ$
Working area length/width	600x600 mm
Standard accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of foot arm rods, a stainless steel sewage basin, a foot controller

BHBD-209 ELECTRIC GYNECOLOGICAL TABLE



Various actions are controlled by Imported LINAK motor

Hand-held remote control and pedal remote control, independent of each other, easy to operate it

With concealed stainless steel sewage waste basin

The mattress is wrapped with imported artificial leather. High density memory sponge is formed in one piece, with no seam

Equipped with brake device which is easy for moving the table

The base is made up of high strength ABS, which is corrosion resistant and easy to clean

Leg plate can be lifted electrically

Concealed stainless steel auxiliary table

With function of auto restoration

SPECIFICATIONS

Model	BHBD-209
Table Length	1270+50 mm
Table width	550+30 mm
Height Adjustment	(700-1000)+50 mm

Trendelenburg/rev-trend	$\geq 5^\circ/\geq 12^\circ$
Backboard downward	$\geq 40^\circ$
leg board reflex	$\geq 30^\circ$
Working area	(350x480)+20 mm
Standard accessory	A pair of arc shaped handle, a pair of leg support, a stainless steel sewage basin, a power cord, a hand held remote controller, a foot controller
Main Voltage AC	220V+22V, 50Hz+1Hz

BHBD-210 ELECTRIC GYNECOLOGICAL TABLE



Adopt high density memory sponge and is once formed, with no seam
 Hydraulic lifting by foot pedal, other actions are controlled by pneumatic
 Imported Y shaped sealing ring and durable
 Stainless steel concealed auxiliary table
 The base is made up of high strength ABS
 Equipped with foot brake device which is easy for moving the table

SPECIFICATIONS

Model	BHBD-210
Table Length	1300+50 mm
Table width	550+30 mm
Height Adjustment	(700 950)+50 mm
Trendelenburg/rev-trend	$\geq 3^\circ/\geq 10^\circ$
Dimension of auxiliary platform	(350x450)+20 mm
Backboard downward	$\geq 70^\circ$
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin

BHBD-211 ELECTRIC GYNECOLOGICAL TABLE



Whole stainless steel hidden auxiliary table with flexible extension and retraction
 Equipped with hidden sewage basin
 With foot control, can adjust the movements easily
 Imported LINAK gear motor, without noise, have a long service life

SPECIFICATIONS

Model	BHBD-211
Table Length	1450+50 mm
Table width	520+30 mm
Height Adjustment	(400 900)+50 mm
Trendelenburg/rev-trend	≥15°/≥25°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Power Supply	220V 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com