



ORTHOPEDIC TRACTION BED

ORTHOPEDIC TRACTION BED

Orthopedic Traction Bed is designed for diagnosis, orthopedic surgeries and traumatic treatment. Our Orthopedic Traction Beds are impeccable on counts of quality, durability and performance, encompasses backrest and leg rest adjustment facilities. Multiple function and the accessories perform expanding traction to different parts of the vertebra and the extremities without moving the human body. It's widely used because of it's easy and convenient operation.

Used in Operation of leg, ICU.

Also known as Orthopedic Bed, Medical bed, Hospital bed.

BHBD-606 ORTHOPEDIC TRACTION BED



SPECIFICATIONS

| | |
|-------------------------|---------------------|
| Model | BHBD-606 |
| Back-rest lifting angle | 0°-70°±5° |
| Foot-rest lifting angle | 0°-30°±5° |
| Load Bearing | 240 kg |
| Height adjustment | 500 mm (from floor) |
| Dimension | 2020x900 mm |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com