



## STAINLESS STEEL BED

## STAINLESS STEEL BED

The Stainless Steel Bed applies to hospital, for nursing, patient rehabilitation treatment and care. These are manufactured using premium quality components in compliance with industrial quality standards. Our hospital beds are based on fundamental principals; safety, ergonomic design, comfort, and easy to use for the patients and hospital staff. Used in Nursing, Patient care.

Also known as Medical bed, hospital bed, Hospital cot.

## BHBD-709 HOSPITAL BED STAINLESS STEEL



## SPECIFICATIONS

Model	BHBD-709
Load Bearing	240 kg
Dimension	2030x900 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)