



ORTHOPEDICS TRACTOR RACK

ORTHOPEDICS TRACTOR RACK

Biolab's Orthopedics tracktor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.

Used in Surgery, Leg operation, Medical unit.

Also known as Operating Table.

400 ORTHOPEDICS TRACTOR RACK



Adopt carbon fiber material, the prespective imaging is clear, light and strong

The plug in design is convient and quick to connect with table

Equipped with a traction frame accessory car, which is convient to move and save space

SPECIFICATIONS

Model	BOPT-401	BOPT-402	BOPT-403
Table Length	1420±50 mm	1700±10 mm	1550±50 mm
Table Width	765±50 mm	750±100 mm	450±50 mm
Traction device horizontal rotary	0-180°		
Traction distance	(0-200)+10%		
Traction Device height adjustment	>400 mm	-	>190 mm



BOPT-401



BOPT-402



BOPT-403



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com