



MULTIFUNCTIONAL OPERATION TABLE

MULTIFUNCTIONAL OPERATION TABLE

Multifunctional operation table are suitable for various operation tables conducted by medical unit such as surgeries of chest, abdomen, gynaecology, obstetrics and orthopaedics.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology.

Also known as Multipurpose operation table, Surgery table, Operating table.

BOPT-501 MULTIFUNCTIONAL OPERATION TABLE



Imported sealing ring with long life

With brake device

all of actions are controlled by both sides

SPECIFICATIONS

Model	BOPT-501
Table Length	2100±50 mm
Table Width	500±20 mm
Height Adjustment	(750~950) ±50 mm
Trendelenburg/Rev-Trend	≥20°/≥15°
Head board downwards	≥90°
Head board upwards	≥40°
Backboard downwards	≥18°
Backboard upwards	≥75°
Leg board downwards	≥90°
Tilt (left)/ Tilt (right)	≥15°/≥15°
Waist board Lifting	100±20 mm
Standard Accessories	A Pair of Shoulder/Waist Support, A pair of Arm Support, A pair of leg Plate, A screen

BOPT-502 MULTIFUNCTIONAL OPERATION TABLE



Imported sealing ring with long life
With brake device
all of actions are controlled by handle under head board

SPECIFICATIONS

Model	BOPT-502
Table Length	2050±50 mm
Table Width	500±20 mm
Height Adjustment	(750~950) ±50 mm
Trendelenburg/Rev-Trend	≥20°
Head board downwards	≥90°
Head board upwards	≥30°
Backboard downwards	≥5°
Backboard upwards	≥75°
Leg board outwards	≥90°
Tilt (left)/ Tilt (right)	≥20°
Waist board Lifting	100±20 mm
Standard Accessories	A Pair of Shoulder/Waist Support, A pair of Arm Support, A pair of leg Plate, A screen



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com