



BOTTLE TOP DISPENSER

BOTTLE TOP DISPENSER

Bottle Top Dispenser combined with routine calibration of equipment yield higher quality results. It makes dispensing liquids easier and the volume is adjusted quickly and set reliably. Low error, capable of reducing time and costs in research and experiment and is fully autoclavable.

Used in Chemical Handling, Pollution Testing, Soil Testing, Environmental Laboratories, Pharmaceutical Laboratories, Research and Development, Qualitative and Quantitative Analysis.

Also known as Laboratory Bottle-Top Dispenser.

BPIP-101 BOTTLE TOP DISPENSER



It is suitable for continuously dispensing liquid with a fixed volume in laboratories and is applicable in various fields

Solvent Resistance: general acid and alkali

Material: ABS

SPECIFICATIONS

Model	BPIP-101
Volume	2.5-25ml
Precision	0.10% (25 μ l)
Accuracy	\pm 0.50% (125 μ l)
Gross Weight (kg)	1.7kg
Dimensions in cubic meters	26.5x17x33cm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com