



MULTI CHANNEL PIPETTE

MULTI CHANNEL PIPETTE

New ergonomic technology and innovative design makes it one of the most innovative pipettes in the market. Multi Channel Pipettes are designed to aspirate and dispense liquid volumes while eliminating many of the potential error risks. Its operation remains correct and controlled, volume settings can be set to prevent accidental volume changes. Used in ELISA, PCR, Cell Culture, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology.

Also known as Laboratory Multi Channel Pipette.

300 MULTI CHANNEL PIPETTE



Light weight

Rotating dispensing head for optimal pipetting convenience.

Due to the optimized standard design of the tip cone, it can be compatible with a wide range of hand of tips.

Large central pipetting button and separate ejection function

True one-handed operation for both right-and left-handers

Completely autoclavable at 121°C(20min)

Volume-change protection

4-Position volume display, always clearly visible externally

Short stroke of only 12.5mm to reduce the risk of RSI(Repetitive Strain Injury)

Corrosion-resistant piston and ejector

Color-coded for easy selection of the right tip

Individual shafts with seals can be easily unscrewed for cleaning or replacing-eliminates expense and long outages.

SPECIFICATIONS

Model	BPIP-301	BPIP-302	BPIP-303	BPIP-304
Autoclavable	121°C			
Volume	0.5-10 µl	5-50 µl	10-100 µl	20-200 µl
Volume Step (µl)	2010-05-01	50/25/5	100/50/10	200/100/20
Increment (µl)	0.01	0.05	0.1	0.2
Accuracy*≤±%	1.6/2/8	0.8/1.4/6	0.8/1.4/4	
CV*≤%	1.0/2/6	0.4/0.8/3	0.3/0.6/2	0.3/0.6/1.5
Type of tips (µl)	20	200	200/300	
Unit	20			

Model	BPIP-305	BPIP-306	BPIP-307	BPIP-308
Autoclavable	121°C	Ejector 121°C		
Volume	30-300 µl	0.5-10 µl	5-50 µl	10-100 µl
Volume Step (µl)	300/150/30	2010-05-01	50/25/5	100/50/10
Increment (µl)	0.3	0.01	0.05	0.1
Accuracy* $\pm\%$	0.6/1.2/3	1.6/2/8	0.8/1.4/6	0.8/1.4/4
CV* $\leq\%$	0.3/0.6/1.5	1.0/2/6	0.4/0.8/3	0.3/0.6/2
Type of tips (µl)	300	20	200	200/300
Unit	20			

Model	BPIP-309	BPIP-310	BPIP-311	BPIP-312
-------	----------	----------	----------	----------

Autoclavable	Ejector 121°C		121°C	
Volume	20-200 µl	30-300 µl	0.5-10 µl	5-50 µl
Volume Step (µl)	200/100/20	300/150/30	2010-05-01	50/25/5
Increment (µl)	0.2	0.3	0.01	0.05
Accuracy* $\pm\%$	0.8/1.4/4	0.6/1.2/3	1.6/2/8	0.8/1.4/6
CV* $\leq\%$	0.3/0.6/1.5		1.0/2/6	0.4/0.8/3
Type of tips (µl)	200/300	300	20	200
Unit	20		5	

Model	BPIP-313	BPIP-314	BPIP-315	BPIP-316
Autoclavable	121°C			Ejector 121°C
Volume	10-100 µl	20-200 µl	30-300 µl	0.5-10 µl
Volume Step (µl)	100/50/10	200/100/20	300/150/30	2010-05-01
Increment (µl)	0.1	0.2	0.3	0.01
Accuracy* $\pm\%$	0.8/1.4/4		0.6/1.2/3	1.6/2/8
CV* $\leq\%$	0.3/0.6/2	0.3/0.6/1.5		1.0/2/6
Type of tips (µl)	200/300		300	20
Unit	5			

Model	BPIP-317	BPIP-318	BPIP-319	BPIP-320
Autoclavable	Ejector 121°C			
Volume	5-50 µl	10-100 µl	20-200 µl	30-300 µl
Volume Step (µl)	50/25/5	100/50/10	200/100/20	300/150/30
Increment (µl)	0.05	0.1	0.2	0.3
Accuracy*≤±%	0.8/1.4/6	0.8/1.4/4		0.6/1.2/3
CV*≤%	0.4/0.8/3	0.3/0.6/2	0.3/0.6/1.5	
Type of tips (µl)	200	200/300		300
Unit	5			

Model	BPIP-321	BPIP-322	BPIP-323	BPIP-324
Autoclavable	121°C			
Volume	0.5-10 µl	5-50 µl	10-100 µl	20-200 µl
Volume Step (µl)	2010-05-01	50/25/5	100/50/10	200/100/20
Increment (µl)	0.01	0.05	0.1	0.2
Accuracy*≤±%	1.6/2/8	0.8/1.4/6	0.8/1.4/4	
CV*≤%	1.0/2/6	0.4/0.8/3	0.3/0.6/2	0.3/0.6/1.5
Type of tips (µl)	20	200	200/300	
Unit	5			

Model	BPIP-325	BPIP-326	BPIP-327	BPIP-328
Autoclavable	121°C			
Volume	30-300 µl	0.5-10 µl	5-50 µl	10-100 µl
Volume Step (µl)	300/150/30	2010-05-01	50/25/5	100/50/10
Increment (µl)	0.3	0.01	0.05	0.1
Accuracy* $\pm\%$	0.6/1.2/3	1.6/2/8	0.8/1.4/6	0.8/1.4/4
CV* $\leq\%$	0.3/0.6/1.5	1.0/2/6	0.4/0.8/3	0.3/0.6/2
Type of tips (µl)	300	20	200	200/300
Unit	5			

Model	BPIP-329	BPIP-330
Autoclavable	121°C	
Volume	20-200 µl	30-300 µl
Volume Step (µl)	200/100/20	300/150/30
Increment (µl)	0.2	0.3
Accuracy* $\pm$ %	0.8/1.4/4	0.6/1.2/3
CV* $\leq$ %	0.3/0.6/1.5	
Type of tips (µl)	200/300	300
Unit	5	



BPIP-301



BPIP-302



BPIP-303



BPIP-304



BPIP-305



BPIP-306



BPIP-307



BPIP-308



BPIP-309



BPIP-310



BPIP-311



BPIP-312



BPIP-313



BPIP-314



BPIP-315



BPIP-316



BPIP-317



BPIP-318



BPIP-319



BPIP-320



BPIP-321



BPIP-322



BPIP-323



BPIP-324



BPIP-325



BPIP-326



BPIP-327



BPIP-328



BPIP-329



BPIP-330



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com