

SEROLOGICAL PIPETTE


SEROLOGICAL PIPETTE

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

600 SEROLOGICAL PIPETTE



Easy to use

Medical grade original Polystyrene

RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation


SPECIFICATIONS

Model	BPIC-601	BPIC-602	BPIC-603	BPIC-604
Volume	1 ml		2 ml	
Precision	2%			
Negative Graduation (ml)	-0.4			
Graduation	1 / 100			
Color	Yellow		Green	
Sterile	Gamma			
Package	Individual	Bulk pack		
Pcs/bag	100 pcs/bag		50 pcs/bag	25 pcs/bag
Pcs/ctn	4000 pcs/carton	3200 pcs/carton	1800 pcs/carton	1200 pcs/carton
Unit	1			


Model	BPIC-605	BPIC-606	BPIC-607	BPIC-608
Volume	5 ml		10 ml	
Precision	2%			
Negative Graduation (ml)	-2.5		-3	
Graduation	10-Jan			
Color	Blue		Red	
Sterile	Gamma			
Package	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag			
Pcs/ctn	800 pcs/carton	600 pcs/carton	4000 pcs/carton	3200 pcs/carton
Unit	1			

Model	BPIC-609	BPIC-610	BPIC-611	BPIC-612
Volume	25 ml		50 ml	
Precision	2%			
Negative Graduation (ml)	-3			
Graduation	05-Jan		1	
Color	Purple		Black	
Sterile	Gamma			
Package	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag			
Pcs/ctn	3000 pcs/carton	2100 pcs/carton	1400 pcs/carton	1000 pcs/carton
Unit	1			


Model	BPIC-613	BPIC-614
Volume	100 ml	
Precision	2%	
Negative Graduation (ml)	-	-
Graduation	1	
Color	Black	
Sterile	Gamma	
Package	Individual	Bulk pack
Pcs/bag	1 pc/bag	
Pcs/ctn	300 pcs/carton	
Unit	1	




BPIC-601




BPIC-602




BPIC-603




BPIC-604




BPIC-605




BPIC-606




BPIC-607




BPIC-608




BPIC-609




BPIC-610




BPIC-611



BPIC-612



BPIC-613



BPIC-614



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com