



TIPS BOX

TIPS BOX

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

700 TIPS BOX



Made from medical grade Polypropylene
Slippy surface with excellent design
DNase-free
Rnase-free
Non-pyrogenic
Filter and non-filter for choice
Bag and rack are available to pack

SPECIFICATIONS

Model	BPIC-701	BPIC-702	BPIC-703	BPIC-704
Volume	10 µl	200 µl	1000 µl	5000 µl
Unit	20			

Model	BPIC-705	BPIC-706	BPIC-707	BPIC-708
Volume	10000 µl	10 µl	200 µl	300 µl
Unit	20			

Model	BPIC-709	BPIC-710	BPIC-711	BPIC-712
Volume	1000 µl	5000 µl	10000 µl	10 µl
Unit	20			10

Model	BPIC-713	BPIC-714	BPIC-715	BPIC-716
Volume	200 µl	300 µl	1000 µl	5000 µl
Unit	10			

Model	BPIC-717	BPIC-718	BPIC-719	BPIC-720
Volume	10000 µl	10 µl	200 µl	300 µl
Unit	10			



BPIC-701



BPIC-702



BPIC-703



BPIC-704



BPIC-705



BPIC-706



BPIC-707



BPIC-708



BPIC-709



BPIC-710



BPIC-711



BPIC-712



BPIC-713



BPIC-714



BPIC-715



BPIC-716



BPIC-717



BPIC-718



BPIC-719



BPIC-720



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com