



UVC STERILIZER

UVC STERILIZER

UVC Sterilizer is cost effective, Chemical free, efficient and environment friendly used to kill bacteria, viruses, mold, fungi and micro-organisms.

Used in Clinics, Disinfection of Microorganisms, Research, Pharmaceuticals.

Also known as Germicidal Lamp.

300 UVC STERILIZER



Simple, easy and convenient-Connected to a wireless network for operation via a mobile phone or tablet devices

Compact and lightweight design-Easy to setup and flexible for any space

Effective disinfection-By using UVC germicidal lamps but no chemicals, it can effectively kill 99.9% germs, also no need to clean up after use

Schedule and timer-Easy to use a mobile phone or tablet to set up timer or schedule for disinfection

SPECIFICATIONS

Model	BSTE-301	BSTE-302	BSTE-303	BSTE-304
Wavelength	253.7 nm			
Lamp Type	UVC Lamp			
Intensity	1360 uW/cm ²	620 uW/cm ²	450 uW/cm ²	400 uW/cm ²
Lamp Holder	2G11 Base			
Wi-Fi	IEEE 802.11 b/g/n 2.4GHz			
Android 4.1, iOS 9.0 or higher version	OS Language			
English, French, Italian, German, Bulgarian, Czech, Danish, Greek, Polish, Portuguese, Russian, Slovak, Spanish, Traditional Chinese, Turkish	Dimension (WxHxD) and Weight(kg)			
260x170x475 mm 2.9 Kg	Package Dimension and Weight	-	-	-
495x235x165 mm 3.3 Kg	Power	-	-	-
144 Watt	Input	-	-	-
260x85x475 mm 1.8 Kg	-	Package Dimension and Weight	-	-
470x209x138 mm 2.4 Kg	-	Power	-	-
72 Watt	-	Input	-	-
260x85x475 mm 1.5 Kg	-	-	Package Dimension and Weight	-
470x209x138 mm 2.2 Kg	-	-	Power	-
36 Watt	-	-	Input	-
260x85x285 mm 1.3 Kg	-	-	-	Package Dimension and Weight

470x209x138 mm 2.0 Kg	-	-	-	Power
18 Watt	-	-	-	Input



BSTE-301



BSTE-302



BSTE-303



BSTE-304



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com