



MOBILE FUME EXTRACTOR

MOBILE FUME EXTRACTOR

The Series of Mobile Fume Extractor is highly efficient purification equipment with self-cleaning system, very ideal for the fume and dust.

BFSD-701 MOBILE FUME EXTRACTOR



Freely moving, easy to operate and install.

The purification system consists of fan, primarily efficient filter, high efficient filter and activated carbon filter.

HEPA filter ensures the filtering efficiency to 99.999% efficiency at 0.3 μ m

SPECIFICATIONS

| | |
|------------------------|---|
| Model | BFSD-701 |
| Air Velocity | 0.3~0.5 m/s |
| Filter Efficiency | 99.995% efficiency at 0.3 μ m |
| Suction Inlet Diameter | 375mm |
| Noise | \leq 75dB(A) |
| External Size (WxDxH) | 500x500x980mm |
| Package Size(WxDxH) | 810x510x600mm(Extracting arm) 630x630x1175mm(Purification system) |
| Gross Weight | 76kg(15kg for Extracting arm; 61kg for Purification system) |
| Net Weight | 67kg |
| Rated Power | 180W |
| Rated Supply | AC110V/220V \pm 10%, 50/60Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com