



## PATHOLOGY WORKSTATION

# PATHOLOGY WORKSTATION

Workstations are used during dissections, autopsies, necropsies, and other pathological procedures. It includes grossing/trimming stations and dissection tables. It may also include a built-in sink and ventilation system. Used in Hospital pathology department, Pathology laboratory. Also known as Pathological sampling bench.

## 300 PATHOLOGY WORKSTATION



- Odor-remove function
- Single water sink, two-person operations
- Exhaust air system: High-power suction device, low noise
- Optional water heater is available, realizing hot & cold water adjustment
- The base cabinet is inward-concave type, convenient for sitting operation
- Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes
- Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning

## SPECIFICATIONS

Model	BAIP-301	BAIP-302	BAIP-303
Work Surface Height	850 mm		
Exhaust Airflow Volume	500~600m3/h		
Main Body Material	304 stainless steel		
Work Table Material	304 stainless steel		
Package Size (WxDxH)	1160x890x2170 mm	1660x890x2170 mm	1960x890x2170 mm
Gross Weight	250 kg	300 kg	340 kg
Working Area Size (WxDxH)	895x430x900 mm	1395x430x900 mm	1695x430x900 mm
External Size (WxDxH)	1000x750x1950 mm	1500x750x1950 mm	1800x750x1950 mm
Consumption	250 W	300 W	
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz		



BAIP-301



BAIP-302



BAIP-303



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)