



UV STERILIZATION CABINET

UV STERILIZATION CABINET

UV Sterilization Cabinet are used to sterile various equipments. They are manufactured to the highest standards to ensure reliable, sterile results.

Used in Hospital, Laboratory,House.

Also known as UV Disinfection Cabinet.

700 UV STERILIZATION CABINET



SPECIFICATIONS

Model	BAIP-701	BAIP-702
Type	With Ozone	Without Ozone
UV Ray Wavelength	253.7/184.9nm	253.7nm
Sterilization Way	2 sides of UV ray+1 side of Ozone	3 sides of UV ray
Capacity	200 L	
Power	120 W	
UV Lamp Material	Quartz glass	
Shelves	Standard 2pc, can be optional	
Internal Size	518x523x778 mm	
External Size	750x663x1200 mm	
Packing Size	875x788x1265 mm	
Gross Weight	100 kg	



BAIP-701



BAIP-702



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com