



## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

## WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## 100 WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



## SPECIFICATIONS

| Model            | BCBT-101                                         | BCBT-102        | BCBT-103        | BCBT-104         |
|------------------|--------------------------------------------------|-----------------|-----------------|------------------|
| Capacity         | 4 gallon (15L)                                   | 12 gallon (45L) | 22 gallon (83L) | 30 gallon (114L) |
| Material         | Double layer anti fire cold-rolled steel PP tray |                 |                 |                  |
| Adjustable Shelf | 1 pc                                             |                 | 3 pcs           | 1 pc             |
| Door             | Single / Manual                                  |                 |                 | Double / Manual  |
| External Size    | 430x430x560 mm                                   | 590x460x890 mm  | 600x460x1650 mm | 1090x460x1120 mm |
| Package Size     | 550x560x725 mm                                   | 710x580x1066 mm | 720x580x1815 mm | 1210x680x1285 mm |
| Gross Weight     | 38 Kg                                            | 66 Kg           | 116 Kg          | 117 Kg           |

| Model            | BCBT-105                                         | BCBT-106         | BCBT-107         |
|------------------|--------------------------------------------------|------------------|------------------|
| Capacity         | 45 gallon (170L)                                 | 60 gallon (227L) | 90 gallon (340L) |
| Material         | Double layer anti fire cold-rolled steel PP tray |                  |                  |
| Adjustable Shelf | 2 pcs                                            |                  |                  |
| Door             | Double / Manual                                  |                  |                  |
| External Size    | 1090x460x1650mm                                  | 860x860x1650 mm  | 1090x860x1650 mm |
| Package Size     | 1210x580x1815 mm                                 | 980x980x1815 mm  | 1210x980x1815 mm |
| Gross Weight     | 155 Kg                                           | 195 Kg           | 243 Kg           |



BCBT-101



BCBT-102



BCBT-103



BCBT-104



BCBT-105



BCBT-106



BCBT-107



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)