



STRONG ACID AND ALKALI STORAGE CABINET

STRONG ACID AND ALKALI STORAGE CABINET

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens.

Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

200 STRONG ACID AND ALKALI STORAGE CABINET



Padlock design, Adjustable shelves, Anti-corrosive materials for strong acid and alkali, Customizable

SPECIFICATIONS

Model	BCBT-201	BCBT-202	BCBT-203	BCBT-204
Capacity	4 gallon (15L)	12 gallon (45L)	20 gallon (75L)	28 gallon (105L)
Material	pp			
Adjustable Shelf	Optional	1 pc		2 pcs
Door	Single / Manual		Double / Manual	
External Size	430x430x500 mm	590x460x890 mm	910x590x890 mm	910x600x900 mm
Package Size	515x515x685 mm	675x545x1025 mm	995x675x1035 mm	995x685x1035 mm
Gross Weight	16 Kg	30 Kg	46 Kg	48 Kg


Model	BCBT-205	BCBT-206	BCBT-207	BCBT-208
Capacity	30 gallon (114L)	45 gallon (170L)	60 gallon (227L)	90 gallon (340L)
Material	pp			
Adjustable Shelf	2 pcs			
Door	Double / Manual			
External Size	1090x460x1120mm	1090x460x1650mm	860x860x1650 mm	1090x860x1650 mm
Package Size	1175x545x1255 mm	1175x545x1785 mm	945x945x1785 mm	1173x943x1785 mm
Gross Weight	58 Kg	72 Kg	90 Kg	107 Kg




BCBT-201



BCBT-202



BCBT-203



BCBT-204



BCBT-205



BCBT-206



BCBT-207



BCBT-208



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com