



CELL COUNTERS

CELL COUNTERS

Bright field detection method with 6 chamber slot per slide
Used in Blood analysis, Diagnostic purposes, Research.
Also known as Cell Counter.

BCLC-301 CELL COUNTER

Bright field detection method with 6 chamber slot per slide

SPECIFICATIONS

Model	BCLC-301
Chamber slot	6 slots per slide
Detection method	Brightfield
Focus	Manual focus
Data transfer	USB
Measuring Time per sample	2s
Sample concentration range	1x 10 ⁴ to 3 x 10 ⁷ cells/ml
Particle/Cell diameter range	5-200 µm
Sample volume	25 µ (*6 slots)
Optics	5 MP CMOS
Display	N/A (connected to laptop)
No. of user profiles	More than 100
Weight	1.8 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com