



DIGITAL DESKTOP ULTRASONIC CLEANER

DIGITAL DESKTOP ULTRASONIC CLEANER

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

100 DIGITAL DESKTOP ULTRASONIC CLEANER



High strength and environmental ABS injection molding
Top quality ultrasonic ceramics parts, quality higher than competitors
Industrial IC, high voltage isolation design, safe using and long working life
Damp proof and anti corrosion process of the PCB
Strict production and quality test control
Digital control, multi stage ultrasonic cleaning adjustment
The tank is polished by SUS304 stainless steel stamping point
2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

SPECIFICATIONS

Model	BULC-101	BULC-102	BULC-103	BULC-104
Capacity	600 ml	750 ml	1400 ml	2500 ml
Working frequency	42 KHz		40 KHz	
Ultrasonic Power	35 W		60 W	70 W
Heating power	No			100 W
Timing time	1-5 mins			1-8 mins
Temperature range	No			65°C
Cover	Yes			
Drain valve	No			
Washing basket	Standard			
Inner tank dimension	155x95x48 mm	150x126x44 mm	168x136x62 mm	247x150x78 mm
Packing size	247x190x145 mm	253x224x150 mm	316x313x213 mm	325x245x210 mm
Gross weight	1.0 Kg	1.17 Kg	2.25 Kg	3.04 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz			



BULC-101



BULC-102



BULC-103



BULC-104



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com