



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

700 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

Model	BULC-701	BULC-702	BULC-703	BULC-704
Capacity	1.3 L	2 L	3 L	6 L
Ultrasonic frequency	40 KHz			
Time range/ Heating temp	1-30 mins / 20- 80°C			
Ultrasonic power	60 W		120 W	180 W
Heating power	200 W			400 W
Tank size	150x140x65 mm	150x140x100 mm	240x135x100 mm	300x150x150 mm
Packing Size	260x260x260 mm	260x260x290 mm	350x260x330 mm	430x280x380 mm
Gross Weight	2.5 Kg	3 Kg	4 Kg	6 Kg

Model	BULC-705	BULC-706	BULC-707	BULC-708
Capacity	10 L	14 L	15 L	19 L
Ultrasonic frequency	40 KHz			
Time range/ Heating temp	1-30 mins / 20- 80°C			
Ultrasonic power	240 W	300 W	360 W	420 W
Heating power	400 W		600 W	
Tank size	300x240x150 mm	300x240x200 mm	330x300x150 mm	330x300x200 mm
Packing Size	430x370x380 mm	430x370x430 mm	460x430x390 mm	460x430x430 mm
Gross Weight	8 Kg	9 Kg	10 Kg	11 Kg

Model	BULC-709	BULC-710	BULC-711	BULC-712
Capacity	22 L	30 L	2 L	
Ultrasonic frequency	40 KHz		80 KHz	
Time range/ Heating temp	1-30 mins / 20- 80°C		1-480 mins / RT-80°C	
Ultrasonic power	480 W	600 W	50 W	

Heating power	800 W	1000 W	100 W
Tank size	500x300x150 mm	500x300x200 mm	150x140x100 mm
Packing Size	630x430x390 mm	630x430x430 mm	260x260x280 mm
Gross Weight	13 Kg	15 Kg	3.5 Kg

Model	BULC-713	BULC-714	BULC-715	BULC-716
Capacity	4 L		6 L	
Ultrasonic frequency	80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	100 W		200 W	
Heating power	400 W			
Tank size	300x150x100 mm		300x150x150 mm	
Packing Size	445x250x400 mm		435x240x430 mm	
Gross Weight	8.5 kg		8.8 Kg	

Model	BULC-717	BULC-718	BULC-719	BULC-720
Capacity	8 L		10 L	
Ultrasonic frequency	80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	200 W		300 W	
Heating power	400 W			
Tank size	300x150x180 mm		300x240x150 mm	
Packing Size	505x355x495 mm		446x346x463 mm	
Gross Weight	13 Kg		12.1 Kg	

Model	BULC-721	BULC-722	BULC-723	BULC-724
Capacity	13 L		27 L	
Ultrasonic frequency	80 Khz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	300 W		500 W	
Heating power	400 W		800 W	
Tank size	300x240x180 mm		500x300x180 mm	
Packing Size	840x670x730 mm		1100x690x820 mm	
Gross Weight	75 kg		99 Kg	

Model	BULC-725	BULC-726	BULC-727	BULC-728
Capacity	27 L			
Ultrasonic frequency	80 Khz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	600 W		700 W	
Heating power	800 W			
Tank size	500x300x180 mm			
Packing Size	1100x690x820 mm			
Gross Weight	101 Kg			



BULC-701



BULC-702



BULC-703



BULC-704



BULC-705



BULC-706



BULC-707



BULC-708



BULC-709



BULC-710



BULC-711



BULC-712



BULC-713



BULC-714



BULC-715



BULC-716



BULC-717



BULC-718



BULC-719



BULC-720



BULC-721



BULC-722



BULC-723



BULC-724



BULC-725



BULC-726



BULC-727



BULC-728



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com