



DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Double adjustable frequencies ultrasonic cleaner uses frequency of 45 KHz /80 KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

800 DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Over- Voltage and over current protection

Standard basket and lid

Low level water protection

Fully microprocessor controller and programmable

Time and temperature real time display

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

Model	BULC-801	BULC-802	BULC-803	BULC-804
Capacity	13 L		27 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	300 W		500 W	
Heating power	400 W		800 W	
Tank size	300x240x180 mm		500x300x180 mm	
Packing Size	840x670x730 mm		1100x690x800 mm	
Gross Weight	73 Kg		101 Kg	

Model	BULC-805	BULC-806	BULC-807	BULC-808
Capacity	27 L			
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins /0- 80°C			
Ultrasonic power	600 W		700 W	
Heating power	800 W			
Tank size	500x300x180 mm			
Packing Size	1100x690x800 mm			
Gross Weight	102 Kg		104 Kg	

Model	BULC-809	BULC-810	BULC-811	BULC-812
Capacity	3 L		4 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	100 W			
Heating power	400 W			
Tank size	230x140x100 mm		230x140x130 mm	
Packing Size	446x250x400 mm		446x346x453 mm	
Gross Weight	7 Kg		12 Kg	

Model	BULC-813	BULC-814	BULC-815	BULC-816
Capacity	6 L		8 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	200 W			
Heating power	400 W			
Tank size	300x150x150 mm		300x150x180 mm	
Packing Size	445x250x430 mm		500x350x480 mm	
Gross Weight	9.5 Kg		15.5 Kg	

Model	BULC-817	BULC-818	BULC-819	BULC-820
Capacity	10 L		13 L	
Adjustable ultrasonic frequency	45 Khz /80 Khz			
Time range/ Heating temp	1-480 mins / RT-80°C			
Ultrasonic power	300 W			
Heating power	400 W			
Tank size	300x240x150 mm		300x240x180 mm	
Packing Size	446x346x483 mm		650x540x630 mm	
Gross Weight	12.5 Kg		37 Kg	



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com