



ULTRASONIC CLEANER

ULTRASONIC CLEANER

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

900 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

Model	BULC-901	BULC-902	BULC-903	BULC-904
Capacity	6.5 L	10 L	15 L	22.5 L
Ultrasonic frequency	Single frequency, 40 KHz			
Time range/ Heating temp	1-99 min/RT-80°C			
Ultrasonic power	180 W	240 W	360 W	480 W
Heating power	500 W			
Tank size	300x150x150 mm	300x240x150 mm	330x300x150 mm	500x300x150 mm
External Size	330x180x310 mm	330x270x310 mm	360x330x310 mm	530x330x310 mm
Packing Size	430x250x350 mm	430x340x350 mm	460x400x400 mm	640x400x400 mm
Gross weight	6 Kg	8 Kg	10 Kg	13 Kg
Power Supply	AC 220V, 50Hz			

Model	BULC-905	BULC-906	BULC-907	BULC-908
Capacity	30 L	6.5 L	10 L	15 L
Ultrasonic frequency	Single frequency, 40 KHz			
Time range/ Heating temp	1-99 min/RT-80°C	1-60 min/RT-80°C		
Ultrasonic power	600 W	180 W	240 W	360 W
Heating power	500 W			
Tank size	500x300x200 mm	300x150x150 mm	300x240x150 mm	330x300x150 mm
External Size	530x330x350 mm	330x180x310 mm	330x270x310 mm	360x330x310 mm
Packing Size	640x400x400 mm	430x250x350 mm	430x340x350 mm	460x400x400 mm
Gross weight	15 Kg	6 Kg	8 Kg	10 Kg
Power Supply	AC 220V, 50Hz			

Model	BULC-909	BULC-910	BULC-911	BULC-912
Capacity	22.5 L	30 L	6.5 L	10 L

Ultrasonic frequency	Single frequency, 40 Khz		Double frequency 28/40 Khz	
Time range/ Heating temp	1-60 min/RT-80°C		1-99 min/RT-80°C	
Ultrasonic power	480 W	600 W	180 W	240 W
Heating power	500 W			
Tank size	500x300x150 mm	300x150x150 mm		300x240x150 mm
External Size	530x330x310 mm	530x330x350 mm	330x180x310 mm	330x270x310 mm
Packing Size	640x400x400 mm		-	-
Gross weight	13 Kg	15 Kg	-	-
Power Supply	AC 220V, 50Hz		-	-

Model	BULC-913	BULC-914	BULC-915	BULC-916
Capacity	15 L	22.5 L	30 L	6.5 L
Ultrasonic frequency	Double frequency 28/40 Khz			Single frequency, 40 Khz
Time range/ Heating temp	1-99 min/RT-80°C			
Ultrasonic power	360 W	480 W	600 W	180 W
Heating power	500 W			
Tank size	330x300x150 mm	500x300x150 mm	500x300x200 mm	300x150x150 mm
External Size	360x330x310 mm	530x330x310 mm	530x330x350 mm	300x180x310 mm
Packing Size	-	-	-	-
Gross weight	-	-	-	-
Power Supply	-	-	-	-

Model	BULC-917	BULC-918	BULC-919	BULC-920
Capacity	10 L	15 L	22.5 L	30 L
Ultrasonic frequency	Single frequency, 40 Khz			
Time range/ Heating temp	1-99 min/RT-80°C			
Ultrasonic power	240 W	360 W	480 W	600 W
Heating power	500 W			
Tank size	300x240x150 mm	330x300x150 mm	500x300x150 mm	500x300x200 mm
External Size	300x270x310 mm	360x330x310 mm	530x330x310 mm	530x330x350 mm
Packing Size	-	-	-	-
Gross weight	-	-	-	-
Power Supply	-	-	-	-



BULC-901



BULC-902



BULC-903



BULC-904



BULC-905



BULC-906



BULC-907



BULC-908



BULC-909



BULC-910



BULC-911



BULC-912



BULC-913



BULC-914



BULC-915



BULC-916



BULC-917



BULC-918



BULC-919



BULC-920



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com