



HORIZONTAL AUTOCLAVE BAHZ-303

HORIZONTAL AUTOCLAVE BAHZ-303

We offer a variety of horizontal laboratory autoclaves for general purpose sterilization applications. Customers to choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

BAHZ-303 HORIZONTAL AUTOCLAVE

- High Quality Stainless steel chamber with drying function
- Safety features for over temperature and pressure
- Equipped with a built in printer
- Automatic opening of safety valve if pressure exceeds 0.24MPa
- Door remains locked if pressure exceeds 0.027MPa
- Over temperature, over pressure, auto protection device
- Automatic power cut off and water cut off in case of water shortage

SPECIFICATIONS

Model	BAHZ-303
Capacity	280 L
Adjustment of Temperature	40°C-134°C
Sterilization Temperature	134°C
Sterilization Time	0-99 min
Heat Average	±2°C
Working Pressure	0.22 MPa
Time of Drying	0-99 min
Chamber Dimension	600x1000 mm
Overall Dimension (LxWxH)	1400x770x1780 mm
Packing Size	960 x 1650 x 2000 mm
Weight	340/440 kg
Power	12000 W
Power Supply	380 V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com