



## UV STERILIZATION CABINET BAIP-701

## UV STERILIZATION CABINET BAIP-701

UV Sterilization Cabinet are used to sterile various equipments.They are manufactured to the highest standards to ensure reliable, sterile results. It rapidly disinfect clinical equipment, furniture and electronic devices.

Used in Hospital, Laboratory,House.

Also known as UV Disinfection Cabinet.

## BAIP-701 UV STERILIZATION CABINET



## SPECIFICATIONS

Model	BAIP-701
Type	With Ozone
UV Ray Wavelength	253.7/184.9nm
Sterilization Way	2 sides of UV ray+1 side of Ozone
Capacity	200 L
Power	120 W
UV Lamp Material	Quartz glass
Shelves	Standard 2pc, can be optional
Internal Size	518x523x778 mm
External Size	750x663x1200 mm
Packing Size	875x788x1265 mm
Gross Weight	100 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)