



**Product Image Coming Soon**

## MEDICAL CLASS N AUTOCLAVE BAMN-401

## MEDICAL CLASS N AUTOCLAVE BAMN-401

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with door locking system, these autoclaves are ideal for small spaces.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## BAMN-401 MEDICAL CLASS N AUTOCLAVE



Product Image Coming Soon

## SPECIFICATIONS

Model	BAMN-401
Capacity	3 L
Sterilization Pressure	0.22 Mp
Sterilization Temperature	134 °C
Sterilization Time	4 min
Sterilization Cycle time	11 min
Required water	Distilled water or purified water
Product Size	300X190X350 mm
Outer Dimension	370x260x460 mm
Gross Weight	7.5 Kg
Power	0.5 KW
Voltage	AC220/110V+10%



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)