



## AUTO COAGULATION ANALYZER BANA-1101

# AUTO COAGULATION ANALYZER BANA-1101

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods  
Used in Clinical diagnostic, Laboratory, Pharmaceutical.  
Also known as coagulation analyzer.

## BANA-1101 AUTO COAGULATION ANALYZER



Optical Coagulation testing theory, PT derived FIB function. With automatic and manual dual mode test functions. Use tilt reagent position, avoid wasting, With LED lighting function. Sampling probe with liquid level sensor and constant heating function, avoid cross contamination

## SPECIFICATIONS

Model	BANA-1101
Analysis Method	Coagulation method
Analysis Items	PT, APPT, TT, FIB
Throughput	140 tests/hour
Sample Positions	5
Tilt reagent position	6
Sample Probe	Liquid level sensor and constant heating function
Calibration	Automatic or manual calibration
Lighting	LED
Interchangeable reagent position	/
Alarm	Abnormal test results alarm, insufficient cleaning fluid alarm, waste liquid overflow alarm
STAT function	Emergency sample available
Software	Windows 7/8/10 or above LIS system available
Power Supply	AC 220V+10%, 60-50Hz; 110V+10% 60Hz, 320W
Package Size	660x640x550 mm
Gross Weight	35 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)