



## AIR SHOWER PASS BOX BAPB-101

## AIR SHOWER PASS BOX BAPB-101

It is necessary passage for people entering clean room in-order to minimize particulate contamination  
Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.  
Also known as Pass Box.

### BAPB-101 AIR SHOWER PASS BOX



Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

## SPECIFICATIONS

Model	BAPB-101
Clean Air Speed	0.36~0.54m/s
UV Lamp	15Wx1
External Size (WxDxH)	660x570x1000 mm
Internal Size (WxDxH)	500x500x500 mm
Material	304 stainless steel
Interlock	Electronical
HEPA Filter	Efficiency ≥99.995% at 0.3μm
Gross Weight	80 kg
Packing Size(WxDxH)	780x700x1250 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)