



AIR SHOWER PASS BOX BAPB-103

AIR SHOWER PASS BOX BAPB-103

It is necessary passage for people entering clean room in-order to minimize particulate contamination
Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.
Also known as Pass Box.

BAPB-103 AIR SHOWER PASS BOX



Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BAPB-103
Clean Air Speed	0.36~0.54m/s
UV Lamp	15Wx1
External Size (WxDxH)	860x770x1200 mm
Internal Size (WxDxH)	700x700x700 mm
Material	304 stainless steel
Interlock	Electronical
HEPA Filter	Efficiency $\geq 99.995\%$ at $0.3\mu\text{m}$
Gross Weight	120 kg
Packing Size(WxDxH)	980x900x1450mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com