



PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-601

PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-601

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

BAPV-601 PRE AND POST VACUUM CLASS B AUTOCLAVE

Class B, Pre & post vacuum
4 preprogrammed cycles, 1- customized programmed cycle
Integrated printer for tractability of sterilization cycles
Automatic regulation of the temperature
Light at the beginning of the cycle; Sound indication of end of cycle
No opening of the door during the cycle; Ability to interrupt the cycle at any time
Malfunction detection and warning system
Safety functions: auto-test program, over heat, over pressure, low water level indicator, safety door lock system present
Standard Accessory Sterilizing tray(2); Waste water tank(1); Tray removal tong(1); Tray holder(1)

SPECIFICATIONS

Model	BAPV-601
Capacity	23 L
Sterilization Temperature	115°C, 121°C, 134°C
Working Pressure	-0.088 Mpa
Vacuum leakage rate	0.13 Kpa/min
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Chamber Dimension	Ø250x450 mm
Overall Dimension (LxWxH)	620x470x385 mm
Display	LCD display
Weight	47 / 51 kg
Power	1900 W
Power Supply	220V/50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com