



PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-801

PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-801

Reliable and robust class B autoclaves intended to deliver highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

BAPV-801 PRE AND POST VACUUM CLASS B AUTOCLAVE

4 preprogrammed sterilization cycles, 1- customized programmed cycles
Built-in Printer & USB port to save and download the recording of sterilization cycle
Open-type water tank for easy cleaning
Light at the beginning of the cycle; Sound indication of end of cycle
No opening of the door during the cycle; Ability to interrupt the cycle at any time
Malfunction detection and warning system, corresponding error code (F1-F19)
Safety: Auto-test program, over heat/over pressure, low water level indicator, safety door lock system ; Safety for draining waste water

SPECIFICATIONS

Model	BAPV-801
Capacity	45 L
Sterilization Temperature	115°C, 121°C, 134°C
Water quality	Distilled water
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Chamber Dimension	Ø320x620 mm
Overall Dimension (LxWxH)	820x585x515 mm
Display	5" color touch screen
Power	3800 W
Power Supply	220V/50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com