



## VERTICAL AUTOCLAVE BAVT-103

## VERTICAL AUTOCLAVE BAVT-103

We offer a variety of vertical laboratory autoclaves for general purpose sterilization applications. Customers can choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Food industry, Chemical industry, Research, Laboratory, Pharmaceutical.

Also known as Vertical Pressure Sterilizer, Vertical Pressure Autoclave, Vertical Sterilizer, Laboratory Vertical Autoclave,.

## BAVT-103 VERTICAL AUTOCLAVE



- Stainless steel structure
- Hand wheel type / Screw type
- Door safety lock system
- Digital Display
- Automatic exhaust system for cold air and pressure
- Over-temperature and over-pressure auto-protection
- Alarm function for completion of process and low water level
- Self-inflating type seal
- Optional: Drying system

## SPECIFICATIONS

Model	BAVT-103
Capacity	75 L
Temperature Range	105~134°C
Heat Average	±1°C
Working Pressure	0.22 MPa
Maximum Working Pressure	0.23 Mpa
Working Temperature	134 °C
Timer Range	0~90 min or 0~99hour59min
Overall Dimension (LxWxH)	550x530x1180 mm
Packing Size	610x610x12580 mm
Weight	113 / 93 kg
Power Supply	3x1.5KW/AC220V 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)