



VERTICAL AUTOCLAVE BAVT-104

VERTICAL AUTOCLAVE BAVT-104

We offer a variety of vertical laboratory autoclaves for general purpose sterilization applications. Customers can choose from the range depending upon their application, volume and budgetary requirements. Features are designed specifically for lab based sterilizations used in research laboratories, pharmaceuticals, microbiology labs, food and chemical industries. High performance, reliable results and a longer life makes it an ideal choice.

Used in Food industry, Chemical industry, Research, Laboratory, Pharmaceutical.

Also known as Vertical Pressure Sterilizer, Vertical Pressure Autoclave, Vertical Sterilizer, Laboratory Vertical Autoclave,.

BAVT-104 VERTICAL AUTOCLAVE



- Stainless steel structure
- Hand wheel type / Screw type
- Door safety lock system
- Digital Display
- Automatic exhaust system for cold air and pressure
- Over-temperature and over-pressure auto-protection
- Alarm function for completion of process and low water level
- Self-inflating type seal
- Optional: Drying system

SPECIFICATIONS

| | |
|---------------------------|---------------------------|
| Model | BAVT-104 |
| Capacity | 100 L |
| Temperature Range | 105~134°C |
| Heat Average | ±1°C |
| Working Pressure | 0.22 MPa |
| Maximum Working Pressure | 0.23 Mpa |
| Working Temperature | 134 °C |
| Timer Range | 0~90 min or 0~99hour59min |
| Overall Dimension (LxWxH) | 580x560x1290 mm |
| Packing Size | 640x640x1360 mm |
| Weight | 125 / 105 kg |
| Power Supply | 3x1.5KW/AC220V 50Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com