



VERTICAL AUTOCLAVE BAVT-3303

VERTICAL AUTOCLAVE BAVT-3303

BAVT-3303 VERTICAL AUTOCLAVE



- Bolt structure
- High quality stainless steel material
- Microcomputer automatic control, arbitrary set of sterilization parameters
- Digital display and touch buttons
- Standard test interface
- The container with automatic protection function: over temperature protection; over pressure protection self discharge; low water level protection, anti dry burning
- With stainless steel bucket or basket
- Buzzer reminder after sterilization, automatic stop.
- Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BAVT-3303
Capacity	75 L
Power	4.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C-134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 400x625
Bucket dimension	Ø 380x560
Basket dimension	Ø 360x280x2
Outer dimension	525x525x1060
Packing size	630x680x1190
Net weight	61 Kg
Gross Weight	82 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com