



VERTICAL AUTOCLAVE BAVT-3306

VERTICAL AUTOCLAVE BAVT-3306

BAVT-3306 VERTICAL AUTOCLAVE



- Bolt structure
- High quality stainless steel material
- Microcomputer automatic control, arbitrary set of sterilization parameters
- Digital display and touch buttons
- Standard test interface
- The container with automatic protection function: over temperature protection; over pressure protection self discharge; low water level protection, anti dry burning
- With stainless steel bucket or basket
- Buzzer reminder after sterilization, automatic stop.
- Automatic cooling air, sterilization end automatic exhaust steam

SPECIFICATIONS

Model	BAVT-3306
Capacity	150 L
Power	6.0 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C-134°C
Sterilization time setting range	4 ~ 120 min
Chamber volume	Ø 500x780
Bucket dimension	Ø 470x720
Basket dimension	Ø 460x350x2
Outer dimension	605x649x1210
Packing size	750x700x1370
Net weight	90 Kg
Gross Weight	116 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com