



## ANALYTICAL BALANCE BBAL-101

## ANALYTICAL BALANCE BBAL-101

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

## BBAL-101 ANALYTICAL BALANCE



Aluminum shell, Anti-jamming, Anti-radiation

The multiple functions modes with tare, counting, unit conversion and weighing

## SPECIFICATIONS

Model	BBAL-101
Capacity (g)	110 g
Resolution (g)	0.001 g
Min Weighing (g)	0.004 g
Stable Time	≤3 sec
Display	LED(Green) / LCD(white back light with black font)
Sensor	Load Cell
Pan Size	Φ 80 / Φ 90 mm
Operation Temperature	5-25°C
Repeatability	±0.002 g
Liner	±0.003 g
Cal. Weight	100 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	7.5 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)