



ANALYTICAL BALANCE BBAL-104

ANALYTICAL BALANCE BBAL-104

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

BBAL-104 ANALYTICAL BALANCE



Aluminum shell, Anti-jamming, Anti-radiation

The multiple functions modes with tare, counting, unit conversion and weighing
It is a magnetic sensor.

SPECIFICATIONS

Model	BBAL-104
Capacity (g)	500 g
Resolution (g)	0.001 g
Min Weighing (g)	0.004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Sensor	Magnet
Pan Size	Φ 80 / Φ 90 mm
Operation Temperature	5-25°C
Repeatability	±0.002 g
Liner	±0.003 g
Cal. Weight	200 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3600606006	RS232C & RS485 Interface
3600606007	Printer



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com